

ROSA[®] Hip System



ROSA[®] Hip

Precisely Personalized

ROSA[®] Hip is a personalized robotic system that enables direct anterior surgeons to evaluate and execute a surgical plan based on real-time feedback and the patient's unique anatomy. It provides surgeons with reassurance and control, while seamlessly integrating into their workflow.

Combining ZBEdge[®] Dynamic Intelligence[™] integrated technologies with the clinical heritage of our comprehensive implant systems, Zimmer Biomet is improving the standard of care.



SURGEON-CENTERED



ACCURATE^{1,2}



EFFICIENT^{3,4}




DATA-DRIVEN

WHY ROBOTICS for Total Hip Arthroplasty?



HIP DISLOCATION IS ONE OF THE LEADING CAUSES OF REVISION⁵⁻⁷ AND CAN INCREASE HOSPITAL COSTS BY UP TO 300% FOLLOWING PRIMARY TOTAL HIP ARTHROPLASTY.^{8,9}



TECHNOLOGIES THAT AID IN COMPONENT IMPLANTATION WILL REDUCE OUTLYING POSITIONS, WHICH MAY REDUCE RISK OF DISLOCATION.¹⁰

*The data refers to increased costs in US hospitals. They may not be applicable to EMEA and vary significantly from country to country.

SURGEON-CENTERED

ROSA Hip is customized to a surgeon's existing direct anterior workflow to maintain surgical efficiency.



FLUOROSCOPY-BASED WORKFLOW
MIRRORS CURRENT PROCEDURE

NO PINS, REFERENCE ARRAYS
OR CAMERAS RESULTING IN NO
ADDITIONAL PROCEDURE SETUP




FLAT PANEL DETECTOR AND IMAGE
INTENSIFIER COMPATIBILITY

Implants Designed to Improve Outcomes


From implants to instruments, surgical systems to support services, each product of the Zimmer Biomet portfolio has been designed to address the distinct needs of individual patients.

The Avenir Complete[®] Hip System, together with the G7[®] Acetabular System, pairs a comprehensive offering of stems, shells and liners (including Dual Mobility) with streamlined instrumentation designed to meet distinct patient needs, simplify the surgical workflow and maximize hospital and OR efficiencies.


ROSA Hip is also compatible with Taperloc[®] Complete Hip System, Avenir[®] Hip System and Echo[®] Bi-Metric[®] Stem.



Avenir Complete[®] Hip System



Taperloc[®] Complete Hip System



Echo[®] Bi-Metric[®] Full Length and Microplasty[®] Hip Stem

ACCURATE^{1,2}

Precisely assists with accurate acetabular component orientation and leg length through robotic guidance.^{1,2}

- Component positioning with ROSA Hip has been shown to be **more accurate and reproducible** than conventional instrumentation^{1,2}
- **Provides real-time data** to evaluate leg length and offset decision-making intra-operatively
 - Designed to measure leg length change from Lesser or Greater Trochanter
- ROSA Hip resulted in **100% of cases within the Lewinnek and Callanan Safe Zones** (fewer outliers compared to conventional instrumentation)¹
- Use **ONE Planner[®] Hip**, our web-based surgical software, to plan a hip replacement case
 - Utilizes pre-operative X-ray images to validate implant components and neck cut location, restore leg length and offset and evaluate pelvic tilt
 - When utilizing ONE Planner Hip, surgeon pre-operative plans will be directly available on the ROSA Hip user interface



EFFICIENT^{3,4}

Simplified setup may minimize additional time to procedure.

- **Auto landmarking and overlay tool supports a streamlined procedure**

- Manual landmarking may be tedious and time consuming. After your first two image uploads, ROSA Hip automatically identifies and places landmarks
- Overlay tool provides a comparison of X-ray images to quickly identify inconsistencies in patient and C-arm positioning
- Ensure femoral rotation is consistent among workflow images with trochanter trace feature

- Plan case virtually with intuitive, web-based **ONE Planner Hip**

- **ONE Matrix™ Hip, our Trial panel**, evaluates the **best possible implant combinations** for each patient in terms of leg length and offset

- No pins, reference arrays or CT scans required

- **Easy to Integrate with Minimal Learning Curve**

- The initial learning curve for the ROSA Hip System can be achieved in 12 cases³
- Robotic-assisted total hip arthroplasty utilizing a fluoroscopy-based system experienced shorter surgical times and OR times compared to CT-based robotic systems⁴



DATA-DRIVEN

Making the best decision when it matters requires data-driven intelligence.

ROSA Hip, a cornerstone of ZBEdge® Dynamic Intelligence™, is an integral part of creating a comprehensive view of orthopedic care informed by data.

Meaningful Connections to Unlock Insights

ZBEdge is Dynamic Intelligence designed to elevate and unlock the full potential of Zimmer Biomet's cutting-edge suite of integrated digital technologies, robotics and implant solutions.



ZBEdge®

by  ZIMMER BIOMET


ZBEdge[®] Analytics

ZBEdge[®] Analytics is a data platform that delivers intra-operative, mobility and outcome insights directly to a smartphone, enabling surgeons to objectively assess their performance and understand the potential impact of clinical decisions on patient recovery.

Within ZBEdge Analytics, ROSA Hip surgeons can monitor their performance and benchmark their progress overtime and versus their peers. Through a personalized dashboard, surgeons can view their individual performance, consistency of surgical times and monitor a set of intra-operative metrics for a cohort of patients.

Data from ZBEdge Analytics is exportable as a CSV. file with one effortless click of a button.

That data is automatically connected, aggregated and segmented across patient populations for ease of use.



ZBEdge[™]
PRIVACY

DATA PRIVACY

**AT ZIMMER BIOMET, THE PATIENT IS ALWAYS THE PATIENT,
AND NEVER THE PRODUCT.**

We accept the responsibility that comes along with the new age of data transformation and we are committed to protecting the patient's privacy.

Our dedicated teams of privacy professionals work to support Zimmer Biomet's data protection obligations, data management and use.

References

1. Kamath AF, Durbhakula SM, Pickering T, Cafferky NL, Murray TG, Wind Jr. MA, Méthot S. Improved Accuracy and Fewer Outliers with a Novel CT-free Robotic THA System in Matchedpair Analysis with Manual THA. *Journal of Robotic Surgery*. 2021 Oct 28. doi: 10.1007/s11701-021-01315-3. Epub ahead of print. PMID: 34709535. Study funded by Zimmer Biomet. Cadaveric studies are not necessarily indicative of clinical performance.
2. Buchan GBJ, Hecht CJ 2nd, Liu D, Mokete L, Kendoff D, Kamath AF. Improved accuracy of a novel fluoroscopy-based robotically assisted THA system compared to manual THA. *J Robot Surg*. 2023 Oct;17(5):2073-2079. doi: 10.1007/s11701-023-01623-w. Epub 2023 May 20. PMID: 37209355.
3. Buchan, G.B.J., et al., The learning curve for a novel, fluoroscopy-based robotic-assisted total hip arthroplasty system. *Int J Med Robot*, 2023; p. e2518.
4. Ong CB, Buchan GBJ, Hecht li CJ, Lawrie CM, DeCook CA, Sculco PK, Kamath AF. Robotic-assisted total hip arthroplasty utilizing a fluoroscopy-guided system resulted in improved intra-operative efficiency relative to a computerized tomography-based platform. *J Robot Surg*. 2023 Dec;17(6):2841-2847. doi: 10.1007/s11701-023-01723-7. Epub 2023 Sep 28. PMID: 37770721.
5. Hailer N P, Weiss R J, Stark A, Kärrholm J. The risk of revision due to dislocation after total hip arthroplasty depends on surgical approach, femoral head size, sex, and primary diagnosis: an analysis of 78,098 operations in the Swedish Hip Arthroplasty Register. *Acta Orthop* 2012; 83 (5): 442-8
6. Kostensalo I, Junnila M, Virolainen P, Remes V, Matilainen3M, Vahlberg T, Pulkkinen P, Eskelinen A, Mäkelä KT. Effect of femoral head size on risk of revision for dislocation after total hip arthroplasty: a population-based analysis of 42,379 primary procedures from the Finnish Arthroplasty Register. *Acta Orthop* 2013; 84 (4): 342-7
7. Gwam CU, Mistry JB, Mohamed NS, Thomas M, Bigart KC, Mont MA, et al. Current epidemiology of revision total hip arthroplasty in the United States: National Inpatient Sample 2009 to 2013. *J Arthroplasty* 2017;32:2088e92
8. Abdel MP, Cross MB, Yasen AT, Haddad FS. The functional and financial impact of isolated and recurrent dislocation after total hip arthroplasty. *Bone Joint J* 2015;97:1046e9.
9. Sanchez-Sotelo J, Haidukewych GJ, Boberg CJ. Hospital cost of dislocation after primary total hip arthroplasty. *J Bone Joint Surg* 2006;88:290.
10. Rowan, Fiachra E., Benjamin, Biju, Pietrak, Jurek R., Haddad, Fares S. Prevention of Dislocation After Total Hip Arthroplasty. *The Journal of Arthroplasty* 33 (2018) 1316e1324.

ROSA

For more information contact
your local representative or visit
zimmerbiomet.com/ROSAHip.



Legal Manufacturer
ROSA Hip
Zimmer CAS
75 Queen Street Suite 3300
Montreal (Quebec) H3C 2N6
Canada
Tel +1.866.3D.ORTHO or
514.395.8883 EC



Legal Manufacturer
ONE Planner Hip
Zimmer CAS
75 Queen Street Suite 3300
Montreal (Quebec) H3C 2N6
Canada
Tel +1.866.3D.ORTHO or
514.395.8883 EC

All content herein is protected by copyright, trademarks and other intellectual property rights, as applicable, owned by or licensed to Zimmer Biomet or its affiliates unless otherwise indicated, and must not be redistributed, duplicated or disclosed, in whole or in part, without the express written consent of Zimmer Biomet.

This material is intended for health care professionals. Distribution to any other recipient is prohibited.

For indications, contraindications, warnings, precautions, potential adverse effects and patient counselling information, see the Instructions for Use or contact your local representative; visit <http://www.zimmerbiomet.eu> for additional product information.

Not for distribution in France. Check for country product clearances and reference product specific instructions for use.

Patients must have a compatible smartphone to use mymobility. Not all patients are candidates for the use of this product and surgeons should evaluate individually to determine which patients are appropriate for remote care. All names used in the mymobility app examples are fictitious. No identification with actual patients or health care professionals is intended or should be inferred. Zimmer Biomet Connected Health mymobility Application. Microsoft and HoloLens are trademarks of the Microsoft Corporation.

© 2024 Zimmer Biomet

 **ZIMMER BIOMET**
Moving You Forward.™

3527.1-EMEA-en-Issue Date-2024-12