
VALO™

ULTRADENT
PRODUCTS, INC.

Improving Oral Health Globally

1. Product Description


With its broadband spectrum, the VALO curing light is designed to polymerize all light cured products in the wavelength range of 385–515 nm per ISO 10650.

The VALO curing light has a medical grade, international power supply and is suitable for power outlets from 100 to 240 volts. The handpiece is designed to rest in a standard dental unit bracket or can be custom mounted using the bracket included with the kit.

Product Components:

- 1 – VALO curing light with 7 ft / 2.1 m cord
- 1 – 9-volt, medical grade, international power supply with 6 ft / 1.8 m cord and international plugs
- 1 – VALO Barrier Sleeve sample pack
- 1 – VALO light shield
- 1 – Curing light surface mounting bracket with double stick adhesive tape

Overview of Controls:

The manufacturer accepts no liability for any damage resulting from the improper use of this unit and/or for any purpose other than those covered by these instructions. For all products described, carefully read and understand all instructions and SDS information prior to use.

2. Indications for Use/Intended Purpose

The source of illumination for curing photo-activated dental restorative materials and adhesives.

3. Warnings and Precautions

Risk Group 2

CAUTION UV emitted from this product. Eye or skin irritation may result from exposure. Use appropriate shielding.

CAUTION Possibly hazardous optical radiation emitted from this product. Do not stare at operating lamp. May be harmful to the eyes.

- DO NOT look directly into the light output. Patient, clinician, and assistants should always wear amber colored UV eye protection when the VALO light is in use.
- To prevent the risk of electric shock, no modification of this equipment is allowed. Use only the included Ultradent VALO power supply and plug adapters. If these components are damaged, do not use and call Ultradent Customer Service to order a replacement.
- Portable RF communications equipment may degrade performance if used closer than 30 cm (12 in.)
- Use only authorized accessories, cables, and power supplies to prevent improper operation, increased electromagnetic emissions, or decreased electromagnetic immunity (refer to Electromagnetic Emissions section)
- To prevent the risk of thermal irritation or injury, avoid back-to-back curing cycles and do not expose oral soft tissues at close proximity for more than 10 seconds in any mode. If longer curing times are required, use multiple curing cycles with periods of rest between cycles or use a dual-cure product to avoid heating soft tissue.
- Use caution when treating patients who suffer from adverse photobiological reactions or sensitivities, patients who are undergoing chemotherapy treatment, or patients being treated with photosensitizing medication
- This unit may be susceptible to strong magnetic or static electric fields, which could disrupt the programming. If you suspect this has occurred, unplug the unit momentarily and then re-plug it into the outlet.
- DO NOT wipe down the VALO curing light with caustic or abrasive cleaners, autoclave, or immerse in any kind of ultrasonic bath, disinfectant, cleaning solution, or liquid. Failure to follow included processing instructions may render device inoperable.
- To help prevent cross contamination and help keep dental composite material from adhering to the surface of the lens and wand body, a barrier sleeve must be used over the VALO light with each use
- To prevent the risk of cross-contamination, barrier sleeves are single patient use
- To reduce the risk of corrosion, remove barrier sleeve after use
- To reduce the risk of under-cured resins, do not use curing light if lens is damaged

4. Stepwise Instructions

Preparation

1. Connect the 9-volt power cord to the handpiece cord.
2. Plug the power cord into any electrical outlet (100-240 VAC). The VALO light handpiece will beep twice when powering on, and the timing lights will illuminate indicating the light is ready for use.
3. Prior to each use, place a new barrier sleeve over the curing light, and minimize wrinkles over the lens for best results.
 - To help prevent cross-contamination and help keep dental composite material from adhering to the surface of the lens and wand body, an Ultradent approved barrier sleeve must be used over the VALO curing light with each use. Barrier sleeves are intended for single-patient use.

VALO Light Shield:

- The VALO light shield is oval-shaped, can be rotated for maximum use, and can be used with a transparent barrier sleeve.

Use

1. Each power mode is used for the curing of dental materials with photo initiators. See Quick Mode Guide for recommended curing times. NOTE: The curing light is programmed to cycle from the Standard Power to the High Power to the Xtra Power mode in sequence. For example, to change from the Standard Power mode to the Xtra Power mode, it is necessary to cycle into the High-Power mode and then to the Xtra Power mode.
2. The curing light stores the most recently used timing interval and mode, and it will default back to this whenever the modes are changed or if the batteries are removed.

Operation

CURING MODE: Standard Power mode

TIMING INTERVALS: 5, 10, 15, and 20 seconds.

- The curing light defaults to this mode when it is INITIALLY powered on. The Mode/Status Light will be green and the four green Timing Lights are illuminated, indicating Standard Power mode.
- To change timing intervals, quickly press the Time/Mode Button.
- Press the Power Button to cure. To stop curing prior to completion of a timing interval, press the Power Button again.

NOTE: 20 sec will deliver the most energy of any mode, See Light Intensity section in the Specification Information Table for energy values.

CURING MODE: High Power mode

TIMING INTERVALS: 1, 2, 3, and 4 seconds.

- From Standard Power mode, press and hold the Time/Mode Change Button for 2 seconds and release. The Mode/Status Light will be orange, and the four green Timing Lights will illuminate and flash, indicating High Power mode.
- To change timing intervals, quickly press the Time/Mode Button.
- Press the Power Button to cure. To stop curing prior to the completion of a timing interval, press the Power Button again.
- To return to Standard Power mode, press and hold the Time/Mode Change Button for 2 seconds and release, this will cycle to Xtra Power mode. Press and hold again for 2 seconds, and release. The Mode/Status Light will be green and the four green Timing Lights are illuminated, indicating Standard mode.

CURING MODE: Xtra Power mode

TIMING INTERVAL: 3 seconds only (Note: The Xtra Power mode has a 2 second safety delay at the end of each curing cycle to limit heating during consecutive curing. At the end of the delay, beeping indicates unit is ready for continued use).

- From Standard Power mode, press the Time/Mode Change Button for 2 seconds, release, press and hold again for 2 seconds, and release. The Mode/Status Light will be orange and flash, and three of the green Timing Lights will illuminate and flash, indicating Xtra Power mode.
- Press the Power Button to cure. To stop curing prior to the completion of a timing interval, press the Power Button again.
- To return to the Standard Power mode, press and hold the Time/Mode Button for 2 seconds and release. The Mode/Status Light will be green and the green Timing Lights are illuminated, indicating Standard Power mode.

Sleep Mode: The curing light will go into SLEEP mode after 1 hour of inactivity, as indicated by a slow flashing of the mode/status light. Pressing any button will wake up the curing light and automatically return it to the last setting used.

Cleanup

1. Discard used barrier sleeves in standard waste after each patient.
2. See Processing Section.

Mounting Bracket Instructions

1. Bracket should be mounted to a flat, oil-free surface.
2. Clean surface with rubbing alcohol.
3. Peel backing off the bracket's adhesive tape.
4. Position bracket so the curing light lifts upward when removed. Press firmly into place.

Quick Mode Guide

Mode	Standard Power	High Power	Xtra Power
------	----------------	------------	------------

Power Button				
Mode/Timing LEDs				
Time Buttons				
Time Options	5s 10s 15s 20s	1s 2s 3s 4s	3s Only	
To Change Time	Press and release Time Button quickly to cycle through time options.			
To Change Modes	Press and hold Time Button for 2 seconds and release. VALO will cycle to next Mode.			
Legend	Solid LEDs		Blinking LEDs	

Quick Curing Guide

Curing Mode	Power* (mW)	Irradiance* (mW/cm2)	Total Exposure Time (Seconds)	Energy** (Joules)
Standard	670	900	5/10/15/20	3.4J/6.7J/10.1J/13.4J
High Power Plus	970	1,300	1/2/3/4	1.0J/1.9J/2.9J/3.9J
Xtra	1,570	2,100	3	4.7J

*Nominal output at a distance of 2 mm from the tip of the lens to the top surface of the composite.

**The values in the total energy per cycle (joules) section are rounded up to the nearest tenth.

Quick Warning Guide:

Warnings	
Call Customer Service for Repair	Call Customer Service for Repair
<ul style="list-style-type: none"> •No sound •Flashing, 2 seconds •Allows operation 	<ul style="list-style-type: none"> •Continuous 3 beeps •Prohibits operations

5. Maintenance

General Cleaning of Curing Light

After each use, moisten a gauze or soft cloth with an approved surface disinfectant and wipe the surface and lens. Unauthorized cleaners may cause damage to the curing light.

ACCEPTABLE CLEANERS:

- 70% Isopropyl alcohol
- 70% Ethanol

General Cleaning of the Light Shield:

Clean the VALO Cordless Light sShield using any surface disinfectant. DO NOT autoclave.

User-Performed Maintenance

1. Use a barrier sleeve to prevent dental composite from adhering to the surface of the lens. If necessary use a plastic or stainless steel dental instrument to carefully remove any adhered composite. Do not use tools that will damage the lens.
2. Light meters differ greatly and are designed for specific light guide tips and lenses. Ultradent recommends routinely checking the output in Standard Power mode. NOTE: the true numeric output will be skewed due to the inaccuracy of common light meters and the custom LED pack in the curing light.

Manufacturer Repair

1. Repairs are only to be performed by authorized service personnel. Ultradent to provide service personnel with documentation to perform repairs.

Warranty

Ultradent Products, Inc. ("Ultradent") warrants that this product shall, for a period of 5 years from the date of purchase, when operated according to the operating instructions included with the product, (i) conform in all material respects to the specifications set forth in Ultradent's documentation accompanying the product; and (ii) be free from defects in material and workmanship.

This limited warranty is nontransferable and applies solely to the original purchaser and does not extend to subsequent owners of the product. This limited warranty does not cover any other accessory components such as, but not limited to, batteries, chargers, adapters, or adaptive lenses. This limited warranty is void if the product fails or is damaged due to negligence, abuse, misuse, accident, modification, tampering, alteration, or failure to follow the applicable instructions for use. For purposes of example only, a product that is dropped and damaged is not covered under this warranty. To qualify under this limited warranty, proof of purchase (e.g., sales receipt or similar documentation) must be submitted to Ultradent along with the defective product.

A defective product meeting the warranty conditions set forth herein will, at Ultradent's sole discretion, either be repaired or replaced. In no event shall Ultradent's liability for the product exceed the purchase price paid by the purchaser. Under no circumstances shall Ultradent be liable for any indirect, incidental, foreseen, unforeseen, special, or consequential damages arising out of or in connection with the use of this product.

6. Processing

- This product is intended to be used with a barrier sleeve. See above for maintenance and general cleaning

7. Storage and Disposal

Curing Light Storage and Transport:

- Temperature: +10°C to +40°C (+50°F to +104°F)
- Relative Humidity: 10% to 95%
- Ambient Pressure: 500 hPa to 1060 hPa

When disposing of electronic waste, (i.e. devices, chargers, batteries and power supplies), follow local waste and recycling guidelines.


8. Technical Considerations

Accessories

Item	CE Information
VALO Light Shield	

Technical information/Data

Effective Composite Curing Wavelength Bands:

Attribute	Information/Specification		
Lens	Diameter 9.75 mm		
Wavelength range	<ul style="list-style-type: none"> Utilizable wavelength range: 385 – 515nm Peak wavelengths: 395 – 415nm and 440 – 480nm 		
Light Intensity Table	Nominal Radiant Exitance Optic Comparison Chart		
	Measurement Instrument	† Gigahertz spectrum analyzer	
		Exitance	Total Power
	Aperture of Meter	15 mm	15 mm
	Standard Power (±10%)	900 mW/cm ²	670 mW
	High Power Plus (±10%)	1300 mW/cm ²	970 mW
Xtra Power (±10%)	2100 mW/cm ²	1570 mW	
VALO Curing Light	Ratings: IEC 60601-1 (Safety), IEC 60601-1-2 (EMC)	Weight: 8 oz./226 g (with cord) Length: 9.26 in./23.5 cm Width: .79 in./2 cm Cord length: 7 ft/2.1 m	
Power Supply	Output - 9VDC at 2A Input - 100VAC to 240VAC Ultrafeed P/N 5930 VALO Power Supply with international Plugs	Rating: IEC 60601-1 (Safety) Cord length - 6 ft/1.8 m VALO power supply provides isolation from MAINS power	

Operating Conditions	Temperature: +10°C to +32°C (+50°F to +90°F) Relative Humidity: 10% to 95% Ambient Pressure: 700 hPa to 1060 hPa
Duty Cycle:	The curing light is designed for short-term operation. At maximum ambient temperature (32°C) 1 minute ON back-to-back cycling, 30 minutes OFF (cooling-off period).

Trouble Shooting

If the solutions suggested below do not rectify the problem, please call Ultradent at 800.552.5512. Outside the United States, call your Ultradent distributor or dental dealer.	
Problem	Possible Solutions
Light will not turn on	<ol style="list-style-type: none"> 1. Press the Time/Mode Change Button or Power Button to wake from Power Save Mode. 2. Check that both cords are firmly connected together and to the electrical outlet. 3. Confirm power to the wall outlet.
Light does not stay on for desired time	<ol style="list-style-type: none"> 1. Check Mode and Timing lights for correct time input. 2. Confirm all cord connections are fully seated. 3. Unplug and re-plug power cord into the electrical receptacle.
Light is not curing resins properly	<ol style="list-style-type: none"> 1. Check lens for residual cured resins/composites. 2. Using proper amber UV eye protection, verify the LED lights are working. 3. Check power level with light meter. If using a light meter, Ultradent recommends checking the VALO light in Standard Power mode. <p>NOTE: The true numeric output will be skewed due to the inaccuracy of common light meters and the custom LED pack the VALO light uses. Light meters differ greatly and are designed for specific light guide tips and lenses.</p> <ol style="list-style-type: none"> 4. Check expiration date on curing resin. 5. Ensure proper technique is being followed (adhesive/composite) by manufacturer recommendations.
Cannot change mode or time intervals	Hold both Time/Mode and Power buttons down until a series of beeps indicates the curing light is unlocked.

9. Miscellaneous Information

Guidance and Manufacturer's Declaration for Electromagnetic Emissions		
The VALO light is intended for use in the electromagnetic environment specified below. The customer or user should ensure that it is used in such an environment. WARNING: Use only authorized accessories, cables, and power supplies to prevent improper operation, increased electromagnetic emissions or decreased electromagnetic immunity.		
Emissions Test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The VALO light uses a Globtek medical grade 9VDC adaptor, operates with brown-out protection, and provides limited EMI, RF, and surge suppression.
RF emission CISPR 11	Class B	The VALO light uses electrical and electromagnetic energy only for their internal functions. However, any RF emissions are very low and are not likely to cause interference in nearby electronic equipment.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/flicker emissions IEC 61000-3-3	COMPLIES	The VALO light is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings for domestic use.

Guidance and Manufacturer's Declaration for Electromagnetic Immunity			
The VALO light is intended for use in the electromagnetic environment specified below. The customer or user should ensure that it is used in such an environment.			
IMMUNITY test	IEC 60601 test level	Compliance level	Electromagnetic environment guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 8 kV contact ± 15 kV air	± 8 kV contact ± 15 kV air	Physical environment should be restricted to the following: 1. IP Code: IP20 2. Do not immerse in liquid. 3. Do not use around flammable gas. Unit is Non-APG and Non-AP. 4. Storage humidity range: 10% - 95% 5.Storage temperatures range: 10° C - 40° C
Electrical fast transient/ burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	± 2 kV for power supply lines Note 1: VALO has no I/O ports	Mains power quality should be that of a typical residential, commercial, or hospital, or military environment
Surge IEC 61000-4-5	± 1 kV line to line ± 2 kV line to earth	± 1 kV line to line ± 2 kV line to earth	
Voltage, dips, shorts, interruptions and variations on the power supply input lines IEC 61000-4-11	<5% U (>95% dip in U for 0.5 cycle) 40% U (60% dip in U for 5 cycles) 70% U (30% dip in U for 25 cycles) <5% U (>95% dip in U for 5 s)	<5% U (>95% dip in U for 0.5 cycle) 40% U (60% dip in U for 5 cycles) 70% U (30% dip in U for 25 cycles) <5% U (>95% dip in U for 5 s) Note 2: Self recovers	Mains power quality should be that of a typical residential, commercial, hospital, or military environment. The Globtek 9VDC medical grade adaptor that is supplied with the VALO light operates from mains ranging from 100VAC – 240VAC and is capable of limited brown out, EMI, and surge protection. If the VALO light user requires continued operations without mains interruption, or the mains in any particular region of a country are considered bad due to continuous brown-out, black-out, or excessively noisy power conditions, it is recommended that the VALO light be powered from an uninterruptible power supply or the customer purchase a VALO™ Cordless curing light unit
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical, residential, home health care, commercial, hospital, or military environment.
NOTE: U is the a.c. mains voltage prior to application of the test level Note 1: The VALO light is not equipped with any ports or any accessible I/O lines. Note 2: If there is a drop of 95% in Mains voltage the VALO light will not operate. It has no internal energy storage mechanism. The VALO light will turn off. When power levels are restored, VALO will restart and return to the same state before power loss. The VALO light will self-recover.			

Guidance and Manufacturer's Declaration for Electromagnetic Immunity for non-life support systems

The VALO light is intended for use in the electromagnetic environment specified below. The customer or user should ensure that it is used in such an environment.

IMMUNITY test	IEC 60601 test level	Compliance level	Electromagnetic environment guidance
Conduction RF	3 Vrms	3 Vrms	Portable and mobile RF communications equipment should be used no closer to any part of the VALO light, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = \left[\frac{3,5}{V_1} \right] \sqrt{P}$ $d = \left[\frac{3,5}{E_1} \right] \sqrt{P} \quad \text{80 MHz to 800 MHz}$ $d = \left[\frac{7}{E_1} \right] \sqrt{P} \quad \text{800 MHz to 2.5 GHz}$ <p>P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range.</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol: </p>
IEC 61000-4-6	150 kHz to 80 MHz	150 kHz to 80 MHz	
Radiated RF	3 V/m	3 V/m	
IEC 61000-4-3	80 MHz to 2.5 GHz	80 MHz to 2.5 GHz	

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people. Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the VALO light is used exceeds the applicable RF compliance level above, the VALO light should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the VALO light.

b. Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Guidance and Manufacture's Declaration for recommended separation distances between portable and mobile RF communications equipment and the VALO curing light

The VALO light is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The user of the VALO light can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the VALO light as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter (P in Watts)	Separation distance according to frequency of transmitter (meters)		
	150 kHz – 80 MHz $d = \left[\frac{3.5}{V_1} \right] \sqrt{P}$	80 MHz – 800 MHz $d = \left[\frac{3.5}{E_1} \right] \sqrt{P}$	800 MHz – 2.5 GHz $d = \left[\frac{7}{E_1} \right] \sqrt{P}$
0.01	0.12 meters	0.035 meters	0.07 meters
0.1	0.37 meters	0.11 meters	0.22 meters
1	1.7 meters	0.35 meters	0.7 meters
10	3.7 meters	1.11 meters	2.22 meters
100	11.7 meters	3.5 meters	7.0 meters

The VALO light has been tested according to IEC 60601-1-2:2014 and passed under radiated field strengths of 10 V/m between 80-MHz to 2.5 GHz. The value of 3Vrms corresponds to V1 and the value 10 V/m corresponds to E1 in the formulas above.

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

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
10%

1060 hPa

MEDICAL - GENERAL MEDICAL EQUIPMENT
AS TO ELECTRICAL SHOCK, FIRE, AND
MECHANICAL HAZARDS ONLY
IN ACCORDANCE WITH
ANSI/AAMI ES 60601-1 (2005) + AMD (2012)
CAN/CSA-C22.2 NO. 60601-1 (2008) +
(2014)
E357088

Rx

EN - For professional use only
DE - Nur zur Anwendung durch den Zahnarzt
FR - Pour usage professionnel seulement
NL - Alleen voor professioneel gebruik
IT - Solo per uso professionale
ES - Solo para uso profesional
PT - Apenas para utilização profissional
SV - Endast för professionell användning
DA - Kun til professionel brug
FI - Vain ammattikäyttöön
EL - Για επαγγελματική χρήση μόνο
CS - Pouze pro odborné použití
BG - Само за професионално ползване
HR - Samo za profesionalnu uporabu
ET - Ainult professionaalseks kasutuseks
HU - Kizárólag professzionális felhasználásra
GA - Le haghaidh úsáide gairmiúla a mhain
LV - Tikai profesionālai lietošanai
LT - Tik profesionaliam naudojimui
MT - Għall-użu professjonali biss
RO - Numai pentru uz profesional
SK - Len na profesionálne použitie
SL - Samo za profesionalno uporabo
NO - Kun for profesjonell bruk
PL - Wyrób wyłącznie do użytku profesjonalnego.
RU - Только для профессионального использования
TR - Sadece profesyonel kullanım için
AR - للاستخدام من قبل المهنيين فقط
TH - สำหรับการใช้งานโดยมืออาชีพเท่านั้น
KO - 전문가 전용
ZH - 仅供专业人士使用

EN - Importer
DE - Importeur
FR - Importateur
NL - Importeur
IT - Importatore
ES - Importador
PT - Importador
SV - Importör
DA - Importør
FI - Maahantuoja
EL - Εισαγωγέας
CS - Dovozce
BG - Вносител
HR - Uvoznik
ET - Importija
HU - Importőr
GA - Allmhairéoir
LV - Importētājs
LT - Importuotojas
MT - Importatur
RO - Importator
SK - Dovozca
SL - Uvoznik
NO - Importør
PL - Importer
RU - Импортёр
TR - İthalatçı
AR - المستورد
TH - ผู้จำหน่าย
KO - 수입자
ZH - 进口商

EN - Medical Device
DE - Medizinprodukt
FR - Dispositif médical
NL - Medisch instrument
IT - Dispositivo medico
ES - Dispositivo medico
PT - Dispositivo médico
SV - Medicinska naprava
DA - Medicinsk anordning
FI - Lääketieteellinen laite
EL - Ιατροτεχνολογικό προϊόν
CS - Lékařské zařízení
BG - Медицинско изделие
HR - Medicinski uređaj
ET - Meditsiiniseade
HU - Orvostechnikai eszköz
GA - Uirlis leighis
LV - Medicīniskā ierīce
LT - Medicinos prietaisai
MT - Apparatu medicu
RO - Dispozitiv medical
SK - Lékařský přístroj
SL - Medicinska naprava
NO - Medisinsk enhet
PL - Wyrób medyczny
RU - Медицинское Изделие
TR - Tıbbi Cihaz
AR - جهاز طبي
TH - อุปกรณ์ทางการแพทย์
KO - 의료 기기
ZH - 醫療器械

VALO™

For product SDS, reordering, and/or complete description of Ultradent's product line, please visit our website www.ultradent.com
or call 1-800-552-5512 (toll free) or 801-572-4200 (outside of U.S.).
Report any serious incident to the manufacturer and the competent authority.

© Copyright 2024 Ultradent Products, Inc.

ULTRADENT
PRODUCTS, INC.

EC REP Ultradent Products GmbH
Am Westhoyer Berg 30
51149 Cologne Germany

Australian Sponsor:
Ultradent Australia Pty Ltd
22/2 Market St. Sydney, NSW 2000 Australia
1800 29 09 29 | www.ultradent.com.au

Manufactured by

 Ultradent Products, Inc.
505 West Ultradent Drive (10200 South)
South Jordan, UT 84095

Manufactured in the USA from globally sourced materials
99409AR15 111924