VALO™

ORTHO CORDED

**ULTRADENT**  
PRODUCTS, INC.

*Improving Oral Health Globally*


## 1. Product Description

With its broadband spectrum, VALO Ortho is designed to polymerize all light cured products in the wavelength range of 385-515nm per ISO 10650. VALO has a medical grade, international power supply and is suitable for power outlets from 100 to 240 volts. The handpiece is designed to rest in a standard dental unit bracket or can be custom mounted using the bracket included with the kit.

### Product Components:

- 1 – VALO Ortho curing light with 7-foot / 2.1-meter cord
- 1 – 9-volt, medical grade, international power supply with 6-foot / 1.8-meter cord and universal plugs
- 1 – VALO Barrier Sleeve sample pack
- 1 – Black Light Lens
- 1 – VALO Light Shield
- 1 – Curing light surface mounting bracket with double stick adhesive tape

### Overview of Control:

The manufacturer accepts no liability for any damage resulting from the improper use of this unit and/or for any purpose other than those covered by these instructions. For all products described, carefully read and understand all instructions and SDS information prior to use.

## 2. Indications for Use/Intended Purpose

The source of illumination for curing photo-activated dental restorative materials and adhesives.

## 3. Warnings and Precautions

### Risk Group 2

CAUTION UV emitted from this product. Eye or skin irritation may result from exposure. Use appropriate shielding.

CAUTION Possibly hazardous optical radiation emitted from this product. Do not stare at operating lamp. May be harmful to the eyes.

- DO NOT look directly into the light output. Patient, clinician, and assistants should always wear amber colored UV eye protection when VALO is in use.
- To prevent the risk of electric shock, no modification of this equipment is allowed. Use only the included Ultradent VALO power supply and plug adapters. If these components are damaged, do not use and call Ultradent Customer Service to order a replacement.
- Portable RF communications equipment may degrade performance if used closer than 30 cm (12 in”).
- Use only authorized accessories, cables, and power supplies to prevent improper operation, increased electromagnetic emissions or decreased electromagnetic immunity (refer to Electromagnetic Emissions section).
- To prevent the risk of thermal irritation or injury, avoid back-to-back curing cycles and do not expose oral soft tissues at close proximity for more than 10 seconds in any mode. If longer curing times are required, use multiple shorter curing cycles or use a dual-cure product to avoid heating soft tissue.
- Use caution when treating patients who suffer from adverse photobiological reactions or sensitivities, patients who are undergoing chemotherapy treatment, or patients being treated with photosensitizing medication.
- This unit may be susceptible to strong magnetic or static electric fields, which could disrupt the programming. If you suspect this has occurred, unplug the unit momentarily and then re-plug it into the outlet.
- DO NOT wipe down the VALO curing light with caustic or abrasive cleaners, autoclave, or immerse in any kind of ultrasonic bath, disinfectant, cleaning solution, or liquid. Failure to follow included processing instructions may render device inoperable.
- To help prevent cross contamination and help keep dental composite material from adhering to the surface of the lens and wand body, a barrier sleeve must be used over the VALO with each use.
  - To prevent the risk of cross-contamination, barrier sleeves are single patient use.
  - To reduce the risk of corrosion, remove barrier sleeve after use.
  - To reduce the risk of under-cured resins, do not use curing light if lens is damaged.
- This unit is programmed for orthodontic use. If using VALO Ortho for any other dental procedures, DO NOT exceed 10 seconds in the standard power mode. For any other dental procedure: to avoid over heating the pulp use two 10 second cures instead of a full 20 second cure.

## 4. Stepwise Instructions

### Preparation

1. Connect the 9-volt power cord to the handpiece cord.
2. Plug the power cord into any electrical outlet (100-240 VAC). The curing light handpiece will beep twice when powering on, and the timing lights will illuminate indicating the light is ready for use.
3. Place the curing light into a standard dental unit mounting bracket or accessory mounting bracket until ready for use.
4. Prior to each use, place a new barrier sleeve over the curing light.

### Installing Hygienic Barrier Sleeves:

The hygienic barrier sleeve is custom fitted to the curing light and keeps the surface of the curing light clean. The barrier sleeve helps prevent cross contamination, helps keep dental composite material from adhering to the surface of the lens and curing light, and prevents discoloration and corrosion from cleaning solutions.

#### Note:

- Using the hygienic barrier sleeve will reduce the light output by 5-10%. Due to the high output power of the curing light, curing has been shown to be substantially equivalent.
- The curing light must be cleaned and sanitized with appropriate cleaning and/or sanitizing agents after each patient. See section titled Processing.

### VALO Light Shield:

The VALO Cordless Light Shield is oval-shaped, can be rotated for maximum use, and can be used with a transparent barrier sleeve  
Use:

1. Each power mode is used for the curing of dental materials with photo initiators. See Quick Mode Guide for recommended curing times.
- NOTE: The curing light is programmed to cycle from the Standard Power to the Xtra Power mode to the Xtra Power Quadrant mode in sequence. For example, to change from the Standard Power mode to the Xtra Power Quadrant mode, it is necessary to cycle into the Xtra Power mode and then to the Xtra Power Quadrant mode.
2. The curing light stores the most recently used timing interval and mode, and it will default back to this whenever the modes are changed or if the batteries are removed.

### Operation

#### **CURING MODE: Standard Power mode**

TIMING INTERVAL: 5, 10, 15, 20 seconds.

- The curing light defaults to this mode when it is INITIALLY powered on. The Mode/Status Light will be green and the four green Timing Lights are illuminated, indicating Standard Power mode.
- To change timing intervals, quickly press the Time/Mode Button.
- Press the Power Button to cure. To stop curing prior to completion of a timing interval, press the Power Button again.

#### **CURING MODE: Xtra Power mode**

TIMING INTERVAL: 1, 2, 3 seconds

- From Standard Power mode, press the Time/Mode Change Button for 2 seconds and release. The Mode/Status Light will be orange and flash, and three of the green Timing Lights will illuminate and flash, indicating Xtra Power mode.
- To change the timing intervals, quickly press the Time/Mode change button.
- Press the Power Button to cure. To stop curing prior to the completion of a timing interval, press the Power Button again.
- To return to the Standard Power mode, press and hold the Time/Mode Button for 2 seconds, release, hold for 2 seconds and release. The Mode/Status Light will be green and the green Timing Lights are illuminated, indicating Standard Power mode.

#### **CURING MODE: Xtra Power Quadrant mode**

The Xtra Power Quadrant modes is designed to allow the clinician to continuously cure one quadrant, or five teeth without needing to remove the VALO Ortho Cordless hand piece from the oral cavity. Each of the five teeth will be exposed for 3 seconds followed by a 2 second rest.

TIMING INTERVAL: 5 automatically timed bursts of 3 seconds each (Note: there is a 2 second safety delay at the end of each curing cycle).

- From Standard Power mode, press the Time/Mode Change Button for 2 seconds, release, press and hold again for 2 seconds, and release. The Mode/Status Light will be orange and flash, and three of the green Timing Lights will illuminate and flash, indicating Xtra Power Quadrant mode.
- Press the Power Button to cure. To stop curing prior to the completion of a timing interval, press the Power Button again.
- To return to the Standard Power mode, press and hold the Time/Mode Button for 2 seconds and release. The Mode/Status Light will be green and the green Timing Lights are illuminated, indicating Standard Power mode.

Sleep Mode: The curing light will go into SLEEP mode after 1 hour of inactivity, as indicated by a slow flashing of the mode/status light. Pressing any button will wake up the curing light and automatically return it to the last setting used.

### Cleanup

1. Discard used barrier sleeves in standard waste after each patient.
2. See Processing Section.

### Mounting Bracket Instructions

1. Bracket should be mounted to a flat, oil-free surface.
2. Clean surface with rubbing alcohol.
3. Peel backing off the bracket's adhesive tape.
4. Position bracket so the curing light lifts upward when removed. Press firmly into place.

### Quick Mode Guide:

Mode	Standard Power	Xtra Power	Xtra Power Quadrant
Power Button			
Mode/Timing LEDs			
Time Buttons			
Time Options	5s 10s 15s 20s	1s 2s 3s	3s Only
To Change Time	Press and release Time Button quickly to cycle through time options.		
To Change Modes	Press and hold Time Button for 2 seconds and release. VALO Ortho will cycle to next Mode.		
Legend	Solid LEDs ● ● ● ●		Blinking LEDs ✨ ✨ ✨

### Quick Curing Guide:

Recommended Curing Times for Optimal Results with VALO			
Mode	Standard Mode	Xtra Power Mode	Xtra Power Quadrant Mode
Metal & Ceramic Brackets (per bracket)	One 10 second cure	Two 3 second cures	Two 3 second cures
Per Layer	One 10 second cure	One 3 second cure	-
Final Cure	Two 20 second cures	Two 3 second cures	-

Note: Exposure settings and times may need to be adjusted due to composite reactivity, shade, distance from the light lens to the composite, and depth of composite layer. It is up to the dental professional to know the requirement for the material they are using to determine the adequate time and settings.

### Quick Warning Guide:

Warnings	
Call Customer Service for Repair	Call Customer Service for Repair
<ul style="list-style-type: none"> <li>No sound</li> <li>Flashing, 2 seconds</li> <li>Allows operation</li> </ul>	<ul style="list-style-type: none"> <li>Continuous 3 beeps</li> <li>Prohibits operations</li> </ul>

## 5. Maintenance

### Repair

#### User-Performed Repair

1. Routinely check the lens for cured dental resins. If necessary, use a non-diamond dental instrument to carefully remove any adhered resin.
2. Light meters differ greatly and are designed for specific light guide tips and lenses. Ultradent recommends routinely checking the output in Standard Power mode. NOTE: the true numeric output will be skewed due to the inaccuracy of common light meters and the custom LED pack in the curing light.

#### Manufacturer Repair

1. Repairs are only to be performed by authorized service personnel. Ultradent to provide service personnel with documentation to perform repairs.

### Warranty

Ultradent hereby warrants that this instrument shall, for a period of 5 years\*, conform in all material respects to the specifications therefore as set forth in Ultradent's documentation accompanying the product and be free from any defects in materials or workmanship. This warranty applies solely to the original purchaser and is not transferable. All defective products are to be returned to Ultradent. There are no user service components of the VALO system. Tampering with VALO will void its warranty. The VALO warranty does not cover customer damage. For example; if a VALO is misused or dropped and the lens breaks, the customer would be responsible to pay for any necessary repairs.

\*With sales receipt indicating the date of sale to the dentist.

## 6. Processing

After each use, moisten a gauze or soft cloth with an approved surface disinfectant and wipe the surface and lens.

### ACCEPTABLE CLEANERS:

- Lysol Brand III Disinfectant Spray (Recommended)
- Isopropyl alcohol
- Ethyl alcohol-based cleaners
- Lysol® Concentrate (alcohol-based only)

### UNACCEPTABLE CLEANERS - DO NOT USE:

- Strong alkali detergent of any type, including hand soaps and dish soaps
- Bleach-based cleaners (e.g. Clorox™, Sterilox™\*)
- Hydrogen Peroxide based cleaners
- Abrasive Cleaners (e.g. Comet Cleanser™\*\*)
- Acetone or hydrocarbon-based cleaners
- MEK (Methyl Ethyl Ketone)
- Birex™
- Glutaraldehyde
- Quaternary Ammonium Chloride salt-based cleaners
- Cavicide1™\* solution or wipes
- Cavicide™\*\* products (Non-Bleach)\*\*

\*Trademark of a company other than Ultradent

\*\* If used, it may fade the color

### CLEANING THE LIGHT SHIELD:

Cold disinfect the VALO Light Shield using any surface disinfectant. DO NOT autoclave.

## 7. Storage and Disposal

### Curing Light Storage and Transport:

- Temperature: +10°C to +40°C (+50°F to +104°F)
- Relative Humidity: 10% to 95%
- Ambient Pressure: 500 hPa to 1060 hPa

When disposing of electronic waste, (i.e. devices, chargers, batteries and power supplies), follow local waste and recycling guidelines.

## 8. Technical Considerations

### Accessories

Item	CE Information		
VALO Barrier Sleeves	  MDSS GmbH Schiffgraben 41 30175 Hanover Germany	Manufactured by: TIDI Products, LLC. 570 Enterprise Drive Neeenah, WI 54956 Made in USA	Distributed by: Ultradent Products Inc 505 West Ultradent Drive (10200 South) South Jordan, UT 84095 USA
VALO Light Shield			

Attribute	Information/Specification					
Lens	Diameter 9.75 mm					
Wavelength range	<ul style="list-style-type: none"> <li>• Utilizable wavelength range: 385 – 515nm</li> <li>• Peak wavelengths: 395 – 415nm and 440 – 480nm</li> </ul>					
Light Intensity Table	Nominal Radiant Exitance Comparison Chart				Radiant Exitance will vary based on instrument capability, measurement method and light placement. † Demetron radiometers and MARC spectrum analyzers should be used as reference only due to having smaller apertures than VALO curing lights. * Demetron radiometers should be used as reference only due to limitations in power and spectral response. ‡ Radiant Exitance conforms to ISO 10650 when measured with a Gigahertz spectrum analyzer.	
	Measurement Instrument	†* Demetron L.E.D. Radiometer	† MARC spectrum analyzer	‡ Gigahertz spectrum analyzer		
	Aperture of Meter	7 mm	3.9 mm	Exitance		Total Power
	Standard Power (±10%)	1200 mW/cm <sup>2</sup>		1100 mW/cm <sup>2</sup>		820 mW
	Xtra Power (±10%)		3200 mW/cm <sup>2</sup> (+/-20%)	2100 mW/cm <sup>2</sup>		1570 mW
	Xtra Power Quadrant (±10%)		3200 mW/cm <sup>2</sup> (+/-20%)	2100 mW/cm <sup>2</sup>		1570 mW
VALO Ortho Curing Light	Ratings: IEC 60601-1 (Safety), IEC 60601-1-2 (EMC)		Weight: 8 ounces/226 grams (with cord) Length: 9.26 inches/23.5 cm Width: .79 inches/ 2 cm Cord length: 7 feet/2.1 meters			
Power Supply	Output - 9VDC at 2A Input - 100VAC to 240VAC Ultradent P/N 5930 VALO Power Supply with Universal Plugs		Rating: IEC 60601-1 (Safety) Cord length - 6 feet/1.8 meters VALO Ortho power supply provides isolation from MAINS power			
Operating Conditions	Temperature: +10°C to +32°C (+50°F to +90°F) Relative Humidity: 10% to 95% Ambient Pressure: 700 hPa to 1060 hPa					
Duty Cycle:	The curing light is designed for short-term operation. At maximum ambient temperature (32°C) 1 minute ON back-to-back cycling, 30 minutes OFF (cooling-off period).					

## Trouble Shooting

If the solutions suggested below do not rectify the problem, please call Ultradent at 800.552.5512. Outside the United States, call your Ultradent distributor or dental dealer.	
Problem	Possible Solutions
Light will not turn on	<ol style="list-style-type: none"> <li>1. Press the Time/Mode Change Button or Power Button to wake from Power Save Mode.</li> <li>2. Check that both cords are firmly connected together and to the electrical outlet.</li> <li>3. Confirm power to the wall outlet.</li> </ol>
Light does not stay on for desired time	<ol style="list-style-type: none"> <li>1. Check Mode and Timing lights for correct time input.</li> <li>2. Confirm all cord connections are fully seated.</li> <li>3. Unplug and re-plug power cord into the electrical receptacle.</li> </ol>
Light is not curing resins properly	<ol style="list-style-type: none"> <li>1. Check lens for residual cured resins/composites.</li> <li>2. Using proper amber UV eye protection, verify the LED lights are working.</li> <li>3. Check power level with light meter. If using a light meter, Ultradent recommends checking the curing light in Standard Power mode. NOTE: The true numeric output will be skewed due to the inaccuracy of common light meters and the custom LED pack the curing light uses. Light meters differ greatly and are designed for specific light guide tips and lenses.</li> <li>4. Check expiration date on curing resin.</li> <li>5. Ensure proper technique is being followed (adhesive/composite) by manufacturer recommendations.</li> </ol>
Cannot change mode or time intervals	Hold both Time/Mode and Power buttons down until a series of beeps indicates the curing light is unlocked.

## 9. Miscellaneous Information

Guidance and Manufacturer's Declaration for Electromagnetic Emissions		
The curing light is intended for use in the electromagnetic environment specified below. The customer or user should ensure that it is used in such an environment. <b>WARNING:</b> Use only authorized accessories, cables, and power supplies to prevent improper operation, increased electromagnetic emissions or decreased electromagnetic immunity.		
Emissions Test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The curing light uses a Globtek medical grade 9VDC adaptor, operates with brown-out protection, and provides limited EMI, RF, and surge suppression.
RF emission CISPR 11	Class B	
Harmonic emissions IEC 61000-3-2	Class A	The curing light uses electrical and electromagnetic energy only for their internal functions. Therefore, any RF emissions are very low and are not likely to cause interference in nearby electronic equipment.
Voltage fluctuations/flicker emissions IEC 61000-3-3	COMPLIES	The curing light is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings for domestic use.

Guidance and Manufacturer's Declaration for Electromagnetic Immunity			
The curing light is intended for use in the electromagnetic environment specified below. The customer or user should ensure that it is used in such an environment.			
IMMUNITY test	IEC 60601 test level	Compliance level	Electromagnetic environment guidance
Electrostatic discharge (ESD)	± 8 kV contact ± 15 kV air	± 8 kV contact ± 15 kV air	Physical environment should be restricted to the following: 1. IP Code: IP20 2. Do not immerse in liquid. 3. Do not use around flammable gas. Unit is Non-APG and Non-AP. 4. Storage humidity range: 10% - 95% 5. Storage temperatures range: 10° C - 40° C
IEC 61000-4-2			
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	± 2 kV for power supply lines Note 1: the curing light has no I/O ports	Mains power quality should be that of a typical residential, commercial, or hospital, or military environment
Surge IEC 61000-4-5	± 1 kV line to line ± 2 kV line to earth	± 1 kV line to line ± 2 kV line to earth	
Voltage, dips, shorts, interruptions and variations on the power supply input lines  IEC 61000-4-11	<5% U (>95% dip in U for 0.5 cycle)  40% U (60% dip in U for 5 cycles)  70% U (30% dip in U for 25 cycles)  <5% U (>95% dip in U for 5 s)	<5% U (>95% dip in U for 0.5 cycle)  40% U (60% dip in U for 5 cycles)  70% U (30% dip in U for 25 cycles)  <5% U (>95% dip in U for 5 s) Note 2: Self recovers	Mains power quality should be that of a typical residential, commercial, hospital, or military environment.  The Globtek 9VDC medical grade adaptor that is supplied with the curing light operates from mains ranging from 100VAC – 240VAC and is capable of limited brown out, EMI, and surge protection.  If the curing light user requires continued operations without mains interruption, or the mains in any particular region of a country are considered bad due to continuous brown-out, black-out, or excessively noisy power conditions, it is recommended that the curing light be powered from an uninterruptible power supply or the customer purchase a VALO Cordless unit.
Power frequency (50/60 Hz) magnetic field  IEC 61000-4-8	30 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical, residential, home health care, commercial, hospital, or military environment.
NOTE: U is the a.c. mains voltage prior to application of the test level Note 1: The curing light is not equipped with any ports or any accessible I/O lines. Note 2: If there is a drop of 95% in Mains voltage the curing light will not operate. It has no internal energy storage mechanism. The curing light will turn off. When power levels are restored, the curing light will restart and return to the same state before power loss. The curing light will self-recover.			

**Guidance and Manufacturer's Declaration for Electromagnetic Immunity for non-life support systems**

The curing light is intended for use in the electromagnetic environment specified below. The customer or user should ensure that it is used in such an environment.

IMMUNITY test	IEC 60601 test level	Compliance level	Electromagnetic environment guidance
Conduction RF	3 Vrms	3 Vrms	Portable and mobile RF communications equipment should be used no closer to any part of the curing light, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = \left[ \frac{3,5}{V_1} \right] \sqrt{P}$ $d = \left[ \frac{3,5}{E_1} \right] \sqrt{P} \quad 80 \text{ MHz to } 800 \text{ MHz}$ $d = \left[ \frac{7}{E_1} \right] \sqrt{P} \quad 800 \text{ MHz to } 2,5 \text{ GHz}$ P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).  Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range.  Interference may occur in the vicinity of equipment marked with the following symbol: 
IEC 61000-4-6	150 kHz to 80 MHz	150 kHz to 80 MHz	
Radiated RF	3 V/m	3 V/m	
IEC 61000-4-3	80 MHz to 2.5 GHz	80 MHz to 2.5 GHz	

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.  
 NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.  
 A field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the curing light is used exceeds the applicable RF compliance level above, the curing light should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the curing light.  
 b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

**Guidance and Manufacturer's Declaration for recommended separation distances between portable and mobile RF communications equipment and the curing light**

The curing light is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The user of the curing light can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the curing light as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter (P in Watts)	Separation distance according to frequency of transmitter (meters)		
	150 kHz – 80 MHz	80 MHz – 800 MHz	800 MHz – 2.5 GHz
	$d = \left[ \frac{3,5}{V_1} \right] \sqrt{P}$	$d = \left[ \frac{3,5}{E_1} \right] \sqrt{P}$	$d = \left[ \frac{7}{E_1} \right] \sqrt{P}$
0.01	0.12 meters	0.035 meters	0.07 meters
0.1	0.37 meters	0.11 meters	0.22 meters
1	1.7 meters	0.35 meters	0.7 meters
10	3.7 meters	1.11 meters	2.22 meters
100	11.7 meters	3.5 meters	7.0 meters

The curing light has been tested according to IEC 60601-1-2:2014 and passed under radiated field strengths of 10 V/m between 80-MHz to 2.5 GHz. The value of 3Vrms corresponds to V1 and the value 10 V/m corresponds to E1 in the formulas above.

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

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© Copyright 2021 Ultradent Products, Inc. Patent pending.
Ultradent Products, Inc.
Ultradent Products GmbH