VALO™

GRAND

CORDLESS BROADBAND LED CURING LIGHT

ULTRADENT
PRODUCTS, INC.

VALO™

GRAND


Product Description

With its broadband spectrum, the VALO Grand curing light is designed to polymerize all light cured products in the wavelength range of 385-515nm per ISO 10650. The VALO Grand curing light uses the Ultradent VALO rechargeable batteries and battery charger. The curing light is designed to rest in a standard dental unit bracket or can be custom mounted using the bracket included with the kit.

VALO Grand Product Components:

- 1 – VALO Grand curing light
- 4 – Ultradent VALO rechargeable batteries
- 1 – Ultradent VALO battery charger with medical grade 12VDC AC power adapter
- 1 – VALO Barrier Sleeve sample pack
- 1 – VALO Grand light shield
- 1 – Curing Light surface mounting bracket with double stick adhesive tape

Overview of Controls:

For all products described, carefully read and understand all instructions and SDS information prior to use.

Indications for Use/Intended Purpose

The source of illumination for curing photo-activated dental restorative materials and adhesives.

Warnings and Precautions

Risk Group 2
CAUTION UV emitted from this product. Eye or skin irritation may result from exposure. Use appropriate shielding.
CAUTION Possibly hazardous optical radiation emitted from this product. Do not stare at operating lamp. May be harmful to the eyes.

- DO NOT look directly into the light output. Patient, clinician, and assistants should always wear amber colored UV eye protection when curing light is in use.
- To prevent the risk of electric shock, no modification of this equipment is allowed. Use only the included Ultradent VALO light power supply and plug adapters. If these components are damaged, do not use and call Ultradent Customer Service to order a replacement.
- Portable RF communications equipment may degrade performance if used closer than 30 cm (12 in”).
- Use only authorized accessories, cables, and power supplies to prevent improper operation, increased electromagnetic emissions or decreased electromagnetic immunity (refer to Electromagnetic Emissions section).
- To avoid the risk of electrical fire associated with handling of batteries:
 - DO NOT mix rechargeable batteries with non-rechargeable batteries or other battery types.
 - DO NOT attempt to charge non-rechargeable batteries. Charging non-rechargeable batteries can result in injury or property damage.
 - DO NOT autovape or spray batteries, as well battery contacts, charger, or AC power adapter with liquid of any kind. If corrosion appears on the contacts of the battery charger, call Ultradent Customer Service to order a replacement.
 - DO NOT charge batteries around flammable materials.
 - DO NOT keep charger in clinical operation.
- To avoid the risk of injury, DO NOT use batteries that are corroded (rust), dented, emit an odor or fluids, have a torn or missing wrapping, or are otherwise damaged. Call Ultradent Customer Service to order replacement batteries.
- Only use recommended batteries. Alternative batteries can cause malfunction.
- To prevent the risk of thermal irritation or injury, avoid back-to-back curing cycles and do not expose oral soft tissues at close proximity for more than 10 seconds in any mode. If longer curing times are required, use multiple curing cycles with periods of rest between cycles, or use a dual-cure product to avoid heating soft tissue.
- Use caution when treating patients who suffer from adverse photobiological reactions or sensitivities, patients who are undergoing chemotherapy treatment, or patients being treated with photosensitizing medication.
- This unit may be susceptible to strong magnetic or static electric fields, which could disrupt the programming. If you suspect this has occurred, unplug the unit momentarily and then re-plug it into the outlet.
- DO NOT wipe down the curing light with caustic or abrasive cleaners, autovape, or immerse in any kind of ultrasonic bath, disinfectant, cleaning solution, or liquid. Failure to follow included processing instructions may render curing light inoperable.
- To avoid damaging the equipment, DO NOT insert fingers, instruments, or other objects into the battery compartment of the curing light.

- To avoid damaging the equipment, DO NOT attempt to clean the gold contacts, or any part of the battery compartment. Call Ultradent Customer Service if there is a concern.
- To prevent the risk of cross-contamination, barrier sleeves are single use.
- To reduce the risk of corrosion, remove barrier sleeve after use.
- To reduce the risk of under-cured resins, do not use curing light if lens is damaged.

Stepwise Instructions

Preparation

- 1) Charge batteries before using the curing light (See section Battery Maintenance).
 - 2) Prior to each use, place a new barrier sleeve over the curing light, minimizing wrinkles over the lens for best results.
- Note:
- To help prevent cross-contamination and help keep dental composite material from adhering to the surface of the lens and wand body, an Ultradent approved barrier sleeve must be used over the VALO curing light with each use. Barrier sleeves are intended for single-patient use.

VALO Cordless Light Shield:

The VALO Cordless Light Shield is oval-shaped, can be rotated for maximum use, and can be used with a transparent barrier sleeve.

Use

- 1) Each power mode is used for the curing of dental materials with photo initiators. See Quick Mode Guide for recommended curing times.
- NOTE:** The curing light is programmed to cycle from the Standard Power to the High-Power Plus to the Xtra Power mode in sequence. For example, to change from the Standard Power mode to the Xtra Power mode, it is necessary to cycle into the High-Power Plus mode and then to the Xtra Power mode.
- 2) The curing light stores the most recently used timing interval and mode, and it will default back to this whenever the modes are changed or if the batteries are removed.

Operation

CURING MODE: Standard Power Mode

TIMING INTERVALS: 5, 10, 15, and 20 seconds.

- The curing light defaults to this mode when it is INITIALLY powered on. The Mode/Status Light will be green and the four green Timing Lights are illuminated, indicating Standard Power mode.
- To change timing intervals, quickly press the Time/Mode Button.
- Press the Power Button to cure. To stop curing prior to completion of a timing interval, press the Power Button again.

CURING MODE: High Power Plus Mode

TIMING INTERVALS: 1, 2, 3, and 4 seconds.

- From Standard Power mode, press and hold the Time/Mode Change Button for 2 seconds and release. The Mode/Status Light will be orange, and the four green Timing Lights will illuminate and flash, indicating High Power mode.
- To change timing intervals, quickly press the Time/Mode Button.
- Press either Power Button to cure. To stop curing prior to the completion of a timing interval, press either Power Button again.
- To return to Standard Power mode, press and hold the Time/Mode Change Button for 2 seconds and release, this will cycle to Xtra Power mode. Press and hold again for 2 seconds, and release. The Mode/Status Light will be green and the four green Timing Lights are illuminated, indicating Standard mode.

CURING MODE: Xtra Power Mode

TIMING INTERVAL: 3 seconds only (Note: The Xtra Power mode has a 2 second safety delay at the end of each curing cycle to limit heating during consecutive curing. At the end of the delay, beeping indicates unit is ready for continued use).

- From Standard Power mode, press the Time/Mode Change Button for 2 seconds, release, press and hold again for 2 seconds, and release. The Mode/Status Light will be orange and flash, and three of the green Timing Lights will illuminate and flash, indicating Xtra Power mode.
- Press the Power Button to cure. To stop curing prior to the completion of a timing interval, press the Power Button again.
- To return to the Standard Power mode, press and hold the Time/Mode Button for 2 seconds and release. The Mode/Status Light will be green and the green Timing Lights are illuminated, indicating Standard Power mode.

Sleep Mode: The curing light will go into Sleep Mode after 60 seconds of inactivity, as indicated by a slow flashing of the mode/status light. Picking up or touching the unit will wake-up the curing light and automatically return it to the last setting used. To maximize battery life, leave the curing light undisturbed when not in use.

Cleanup

- 1) Discard used barrier sleeves in standard waste after each patient.
- 2) See Processing Section.

Mounting Bracket Instructions

- 1) Bracket should be mounted to a flat, oil-free surface.
- 2) Clean surface with rubbing alcohol.
- 3) Peel backing off the bracket's adhesive tape.
- 4) Position bracket so the curing light lifts upward when removed. Press firmly into place.

Quick Mode Guide:

Mode	Standard Power 1000 mW/cm ²				High Power Plus 1600 mW/cm ²				Xtra Power 3200 mW/cm ²
Power Button									
Mode/Timing LEDs									
Time Buttons									
Time Options	5s	10s	15s	20s	1s	2s	3s	4s	3s Only
To Change Time	Press and release Time Button quickly to cycle through time options.								
To Change Modes	Press and hold Time Button for 2 seconds and release. Curing light will cycle to next Mode.								
Legend	Solid LEDs						Blinking LEDs		

Quick Curing Guide:

Curing Mode	Power* (mW)	Irradiance* (mW/cm ²)	Total Exposure Time (Seconds)	Energy** (Joules)
Standard	970	900	5/10/15/20	4.8J/9.7J/14.5J/19.4J
High Power Plus	1,620	1,500	1/2/3/4	1.6J/3.2J/4.8J/6.5J
Xtra	2,260	2,100	3	6.8J

*Nominal output at a distance of 2 mm from the tip of the lens to the top surface of the composite.

**The values in the total energy per cycle (joules) section are rounded up to the nearest tenth.

Quick Warning Guide:

Power Level Warning	Temperature Warning	Service Warning	LED Warning
Replace batteries	Allow for cool down	Call Customer Service for Repair	Call Customer Service for Repair
<ul style="list-style-type: none"> Low battery: slow flashing Shut off: 3 beeps, flashing Prohibits operation 	<ul style="list-style-type: none"> 3 beeps Slow flashing Prohibits operation 	<ul style="list-style-type: none"> No sound Flashing, 2 seconds Allows operation 	<ul style="list-style-type: none"> Continuous 3 beeps Single rapid flash Prohibits operations

Maintenance

General Cleaning of Curing Light

After each use, moisten a gauze or soft cloth with an approved surface disinfectant and wipe the surface and lens. Unauthorized cleaners may cause damage to the curing light.

ACCEPTABLE CLEANERS:

- 70% Isopropyl alcohol
- 70% Ethanol

General Cleaning of the Light Shield:

Clean the VALO Cordless light shield using any surface disinfectant. DO NOT autoclave.

Repair

User-Performed Maintenance

1) Use a barrier sleeve to prevent dental composite from adhering to the surface of the lens. If necessary, use a plastic or stainless steel dental instrument to carefully remove any adhered composite. Do not use tools that will damage the lens.

2) Light meters differ greatly and are designed for specific light guide tips and lenses. Ultradent recommends routinely checking the output in Standard Power mode. NOTE: the true numeric output will be skewed due to the inaccuracy of common light meters and the custom LED pack in the curing light

Battery Maintenance

Charging and Changing Batteries

The curing light comes with 4 rechargeable lithium iron phosphate batteries. Ultradent recommends using VALO rechargeable batteries only.

How to charge batteries:

- 1) Place charger and batteries in an area where they will not come into contact with liquids or disinfectant spray.
- 2) Insert batteries into charger with positive (+) end pointed towards the indicator lights on charger.
- 3) Green Lights indicate batteries are ready for use.
- 4) Batteries will take 1-3 hours to charge. Leave batteries in charger until ready for use.

NOTE: If the red light on the charger does not turn green when the batteries have been charging longer than three hours, the battery may be degraded (corrupted) and cannot be charged. Try a new battery or call Ultradent Customer Service to order a new set of rechargeable batteries.

How to change/insert batteries:

- 1) Remove back cap by twisting counterclockwise one quarter turn.
- 2) Remove batteries.
- 3) Insert fresh batteries positive (+) side first.
- 4) Reattach back cap by aligning and gently pushing while twisting clockwise. The cap will click when fully attached.
- 5) The unit is ready for use.

Low Batteries: The curing light signals the user that it is time to change the batteries when the Low Battery Indicator Light is flashing red. If the battery charge becomes too low, an audible 3 beep warning sound will occur and the curing light will not allow further operation until batteries are recharged or new batteries are inserted. (See Quick Warning Guide)

Charging Time and Battery Life: Fully charged battery life in the curing light is dependent on the Mode/Time interval, battery type, amount of use, and LED efficiency.

- Recommended Recharge Interval: Recharge as necessary or, when low battery indicator comes on. Life expectancy depends on use.
- Extra Batteries: The curing light comes with 4 rechargeable batteries.

Manufacturer Repair

- 1) Repairs are only to be performed by authorized service personnel. Ultradent is to provide service personnel with documentation to perform repairs.
- 2) When sending units in for repair, service, or calibrations always remove the batteries from the curing light and charger. Wrap batteries, charger, adapter, and the curing light separately in the return box.
- 3) Ship batteries in accordance with local regulations.

Warranty

Ultradent Products, Inc. ("Ultradent") warrants that this product shall, for a period of 5 years from the date of purchase, when operated according to the operating instructions included with the product, (i) conform in all material respects to the specifications set forth in Ultradent's documentation accompanying the product; and (ii) be free from defects in material and workmanship.

This limited warranty is nontransferable and applies solely to the original purchaser and does not extend to subsequent owners of the product. This limited warranty does not cover any other accessory components such as, but not limited to, batteries, chargers, adapters, or adaptive lenses. This limited warranty is void if the product fails or is damaged due to negligence, abuse, misuse, accident, modification, tampering, alteration, or failure to follow the applicable instructions for use. For purposes of example only, a product that is dropped and damaged is not covered under this warranty. To qualify under this limited warranty, proof of purchase (e.g., sales receipt or similar documentation) must be submitted to Ultradent along with the defective product.

A defective product meeting the warranty conditions set forth herein will, at Ultradent's sole discretion, either be repaired or replaced. In no event shall Ultradent's liability for the product exceed the purchase price paid by the purchaser. Under no circumstances shall Ultradent be liable for any indirect, incidental, foreseen, unforeseen, special, or consequential damages arising out of or in connection with the use of this product.

Processing

This product is intended to be used with a barrier sleeve. See above for maintenance and general cleaning.

Storage and Disposal

If storing the curing light for periods longer than 2 weeks, or packing it for travel, always remove the batteries. If batteries are left in the unit for long periods of time without recharging they may become nonfunctional or un-chargeable. Do not store batteries in temperatures over 40°C (104°F) or in direct sunlight.

Curing Light Storage and Transport:

- Temperature: +10°C to +40°C (+50°F to +104°F)
- Relative Humidity: 10% to 95%
- Ambient Pressure: 500 hPa to 1060 hPa

When disposing of electronic waste (i.e. curing lights, chargers, batteries and power supplies), follow local waste and recycling guidelines.

Technical Considerations


Accessories

Item	CE Information
VALO Cordless Light Shield	

Item	CE Information
VALO Charger	
VALO Batteries	

Technical information/Data

Effective Composite Curing Wavelength Bands:

Attribute	Information/Specification			
Lens	Diameter 11.7 mm			
Wavelength range	<ul style="list-style-type: none"> Utilizable wavelength range: 385 – 515nm Peak wavelengths: 395 – 415nm and 440 – 480nm 			
Light Intensity Table	Nominal Radiant Exitance Comparison Chart			Radiant Exitance will vary based on instrument capability, measurement method, and light placement. ‡ Radiant Exitance conforms to ISO 10650 when measured with a Gigahertz-Optik spectrum analyzer.
	Measurement Instrument	‡ Gigahertz spectrum analyzer		
		Exitance	Total Power	
	Aperture of Meter	15 mm	15 mm	
	Standard Power (±10%)	900 mW/cm ²	970 mW	
	High Power Plus (±10%)	1500 mW/cm ²	1615 mW	
Xtra Power (±10%)	2100 mW/cm ²	2260 mW		
VALO Grand Curing Light	Ratings: IEC 60601-1 (Safety), IEC 60601-1-2 (EMC)	Weight: <ul style="list-style-type: none"> With batteries: 6 oz. (190 g) Without batteries: 5 oz. (150 g) Dimension: (8 x 1.28 x 1.06) inches, (203 x 32.5 x 27) mm		
Charger Power Supply	Output - 12VDC at 500mA Input - 100VAC to 240VAC Ultradent P/N 5930 VALO Charger Power Supply with international plug inserts	Ratings: IEC 60601-1 (Safety) Cord Length - 6 ft (1.8 m) VALO Charger Power Supply is a Medical Grade Class II power supply and provides isolation from MAINS power		

Attribute	Information/Specification
VALO Charger	<p>VALO 3.6VDC Lithium Iron Phosphate smart battery charger:</p> <ul style="list-style-type: none"> • Automatic shut off when fully charged • Auto-detection of defective batteries • Protections: Overcharge, Short-circuit, reverse polarity <ul style="list-style-type: none"> ◦ Red LED – Charging ◦ Green LED – Empty or Fully Charged ◦ LED off – short circuit • Charging time: 1 – 3 hours <p>Rating: CE, WEEE</p>
VALO Batteries	<p>Rechargeable: Safe chemistry Lithium Iron Phosphate (LiFePO4) RCR123A</p> <ul style="list-style-type: none"> • Working Voltage: 3.2VDC, 400mAh (Watt hour needed) <p>Ratings: UL 62133-2, CE, RoHS, WEEE</p> <p>Cells are rated for >2000 charge/discharge cycles before hitting 70% of its capacity. It is advisable to replace batteries every 3-5 years, but life may vary depending on use. Replace as necessary.</p>
Operating, Charging Conditions	<p>Temperature: +10°C to +32°C (+50°F to +90°F)</p> <p>Relative Humidity: 10% to 95%</p> <p>Ambient Pressure: 700 hPa to 1060 hPa</p>
Duty Cycle	<p>The curing light is designed for short-term operation. At maximum ambient temperature (32°C) 1 minute ON back-to-back cycling, 30 minutes OFF (cooling-off period).</p>

Trouble Shooting

If the solutions suggested below do not rectify the problem, please call Ultradent at 800.552.5512. Outside the United States, call your Ultradent distributor or dental dealer.	
Problem	Possible Solutions
Light will not turn on	<ol style="list-style-type: none"> 1) Press the Time/Mode or Power Button to wake from Power Save Mode. 2) Check the red Low Battery Indicator for battery charge status. 3) Check that fresh batteries are correctly inserted into the unit. 4) If red and yellow Warning LEDs are flashing this means the curing light has reached its internal temperature safety limit. Allow the curing light to cool down for 10 minutes or use a cool moist towel to cool the unit down quickly. 5) If red Warning LED flashes and beeps continuously, call Ultradent Customer Service for repair.
Light does not stay on for desired time	<ol style="list-style-type: none"> 1) Check Mode and Timing lights for correct time input. 2) Check the Low Battery Indicator for battery charge status. 3) Check that fresh batteries are properly inserted into the unit.
Light is not curing resins properly	<ol style="list-style-type: none"> 1) Check lens for residual cured resins/composites. 2) Using proper amber UV eye protection, verify the LED lights are working. 3) Check power level with light meter. If using a light meter, Ultradent recommends checking the curing light in Standard Power mode. <p>NOTE: The true numeric output will be skewed due to the inaccuracy of common light meters and the custom LED pack the curing light uses. Light meters differ greatly and are designed for specific light guide tips and lenses.</p> <ol style="list-style-type: none"> 4) Check expiration date on curing resin. 5) Ensure proper technique is being followed according to manufacturer recommendations.
Batteries will not charge	<ol style="list-style-type: none"> 1) Make sure batteries are inserted into the charger in the correct orientation and allow batteries to charge for 1-3 hours. 2) If red lights on the charger do not change to green, call Ultradent Customer Service to order replacement batteries and/or charger. 3) If neither green nor red lights on the charger are visible, call Ultradent Customer Service to order or replace charger and/or AC adapter.
Charger does not charge batteries	<ol style="list-style-type: none"> 1) Make sure charger is plugged in and AC adapter is plugged into a working power outlet. 2) If green or red lights on the charger are not visible, call Ultradent Customer Service for new charger and/or AC adapter.
Cannot change mode or time intervals	<ol style="list-style-type: none"> 1) Hold both Time/Mode and Power buttons down until a series of beeps indicates the curing light is unlocked.

Miscellaneous Information

Guidance and Manufacturer's Declaration for Electromagnetic Emissions		
The curing light is intended for use in the electromagnetic environment specified below. The customer or user should ensure that it is used in such an environment.		
Emissions Test	Compliance	Electromagnetic environment - guidance

Guidance and Manufacture's Declaration for Electromagnetic Emissions			
RF emissions CISPR 11	Group 1	The curing light uses batteries and is not affected by EMI, RF, or surge suppression.	
RF emission CISPR 11	Class B	The curing light uses electrical and electromagnetic energy only for their internal functions. However, any RF emissions are very low and are not likely to cause interference in nearby electronic equipment.	
Harmonic emissions IEC 61000-3-2	N/A	Harmonic emission and voltage fluctuation testing are not applicable to the curing light because it is battery powered.	
Voltage fluctuations/flicker emissions IEC 61000-3-3	N/A	The curing light is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings for domestic use.	

Guidance and Manufacture's Declaration for Electromagnetic Immunity			
The curing light is intended for use in the electromagnetic environment specified below. The customer or user should ensure that it is used in such an environment.			
IMMUNITY test	IEC 60601 test level	Compliance level	Electromagnetic environment guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 8 kV contact ± 15 kV air	± 8 kV contact ± 15 kV air	Physical environment should be restricted to the following: 1) IP Code: IP20 2) Do not immerse in liquid. 3) Do not use around flammable gas. Unit is Non-APG and Non-AP. 4) Storage humidity range: 10%-95% 5) Storage temperatures range: 10° C-40° C
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	± 2 kV for power supply lines Note 1: The curing light has no I/O ports	The curing light is battery powered and is not capable of connection to AC MAINS power.
Surge IEC 61000-4-5	± 1 kV line to line ± 2 kV line to earth	± 1 kV line to line ± 2 kV line to earth	Because the curing light is battery powered, it is not subject to electric transients, surges, voltage dips, shorts, interruptions, or variations on AC MAINS power.
Voltage, dips, shorts, interruptions, and variations on the power supply input lines IEC 61000-4-11	<5% U (>95% dip in U for 0.5 cycle) 40% U (60% dip in U for 5 cycles) 70% U (30% dip in U for 25 cycles) <5% U (>95% dip in U for 5 s)	<5% U (>95% dip in U for 0.5 cycle) 40% U (60% dip in U for 5 cycles) 70% U (30% dip in U for 25 cycles) <5% U (>95% dip in U for 5 s) Note 2: Self recovers	The accessory battery charger may be subject to the above but it is separate from and not critical to the operation of the VALO Cordless curing light. If battery voltage of the curing light drops to 4VDC the unit will not allow operation. The VALO Cordless curing light will turn off. When new batteries are put in and the proper power levels are restored, the curing light will restart and return to the same state before power loss. The curing light will self-recover in the event of power loss.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical, residential, home health care, commercial, hospital, or military environment.
Note 1: The curing light is not equipped with any ports or any accessible I/O lines.			

Guidance and Manufacture's Declaration for Electromagnetic Immunity for non-life support systems			
The curing light is intended for use in the electromagnetic environment specified below. The customer or user should ensure that it is used in such an environment.			
IMMUNITY test	IEC 60601 test level	Compliance level	Electromagnetic environment guidance

Guidance and Manufacturer's Declaration for Electromagnetic Immunity for non-life support systems			
Conduction RF	3 Vrms	3 Vrms	Portable and mobile RF communications equipment should be used no closer to any part of the VALO GRAND, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance: $d = \left[\frac{3,5}{V_1} \right] \sqrt{P}$
IEC 61000-4-6	150 kHz to 80 MHz	150 kHz to 80 MHz	$d = \left[\frac{3,5}{E_1} \right] \sqrt{P}$ 80 MHz to 800 MHz
Radiated RF	3 V/m	3 V/m	$d = \left[\frac{7}{E_1} \right] \sqrt{P}$ 800 MHz to 2.5 GHz P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).
IEC 61000-4-3	80 MHz to 2.5 GHz	80 MHz to 2.5 GHz	Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey [†] , should be less than the compliance level in each frequency range [‡] . Interference may occur in the vicinity of equipment marked with the following symbol: 

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.
NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.
[†] Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the curing light is used exceeds the applicable RF compliance level above, the VALO Cordless curing light should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the curing light.
[‡] Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Guidance and Manufacturer's Declaration for recommended separation distances between portable and mobile RF communications equipment and the VALO Cordless curing light			
The curing light is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The user of the curing light can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the curing light as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter (P in Watts)	Separation distance according to frequency of transmitter (meters)		
	150 kHz – 80 MHz $d = \left[\frac{3,5}{V_1} \right] \sqrt{P}$	80 MHz – 800 MHz $d = \left[\frac{3,5}{E_1} \right] \sqrt{P}$	800 MHz – 2.5 GHz $d = \left[\frac{7}{E_1} \right] \sqrt{P}$
0.01	0.12 meters	0.035 meters	0.07 meters
0.1	0.37 meters	0.11 meters	0.22 meters
1	1.7 meters	0.35 meters	0.7 meters
10	3.7 meters	1.11 meters	2.22 meters
100	11.7 meters	3.5 meters	7.0 meters
The curing light has been tested according to IEC 60601-1-2:2014 and passed under radiated field strengths of 10 V/m between 80-MHz to 2.5 GHz. The value of 3Vrms corresponds to V1 and the value 10 V/m corresponds to E1 in the formulas above. For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer. NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies. NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			

Report any serious incident to the manufacturer and the competent authority.
www.ultradent.com / 1-800-552-5512 / 801-572-4200

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For immediate reorder and/or complete descriptions of Ultradent's product line, refer to Ultradent's catalog or call toll free 1-800-552-5512. Outside U.S. call +1-801-572-4200 or visit www.ultradent.com. Report any serious incident to the manufacturer and the competent authority.

800.552.5512 **ULTRADENT.COM**

Manufactured by Ultradent Products, Inc. 505 West Ultradent Drive
(10200 South), South Jordan, UT 84095.
Made in the USA from globally sourced materials.

EC REP

Ultradent Products GmbH
Am Westthor Berg 30
51149 Cologne Germany

Australian Sponsor:
Ultradent Australia Pty Ltd
22/2 Market St. Sydney, NSW 2000
1800 29 09 29 | www.ultradent.com.au

EQUIPMENT

ULTRADENT
PRODUCTS, INC.