
VALO™

C O R D L E S S

ULTRADENT
PRODUCTS, INC.

Improving Oral Health Globally


1. Product Description

With its broadband spectrum, VALO Cordless is designed to polymerize all light cured products in the wavelength range of 385-515nm per ISO 10650. The VALO Cordless curing light uses Ultradent VALO rechargeable batteries and battery charger. The curing light is designed to rest in a standard dental unit bracket or can be custom mounted using the bracket included with the kit.

VALO Cordless Product Components:

- 1 – VALO Cordless curing light
- 4 – Ultradent VALO rechargeable batteries
- 1 – Ultradent VALO battery charger with medical grade 12VDC AC power adapter
- 1 – VALO Barrier Sleeve sample pack
- 1 – VALO Cordless Light Shield
- 1 – Curing light surface mounting bracket with double stick adhesive tape

Overview of Controls:

For all products described, carefully read and understand all instructions and SDS information prior to use.

2. Indications for Use/Intended Purpose

The source of illumination for curing photo-activated dental restorative materials and adhesives.

3. Warnings and Precautions

Risk Group 2

CAUTION UV emitted from this product. Eye or skin irritation may result from exposure. Use appropriate shielding.

CAUTION Possibly hazardous optical radiation emitted from this product. Do not stare at operating lamp. May be harmful to the eyes.

- DO NOT look directly into the light output. Patient, clinician, and assistants should always wear amber colored UV eye protection when curing light is in use.
- To prevent the risk of electric shock, no modification of this equipment is allowed. Use only the included Ultradent VALO power supply and plug adapters. If these components are damaged, do not use and call Ultradent Customer Service to order a replacement.
- Portable RF communications equipment may degrade performance if used closer than 30 cm (12 in”).
- Use only authorized accessories, cables, and power supplies to prevent improper operation, increased electromagnetic emissions or decreased electromagnetic immunity (refer to Electromagnetic Emissions section).
- To avoid the risk of electrical fire associated with handling of batteries:
 - o DO NOT mix rechargeable batteries with non-rechargeable batteries or other battery types.
 - o DO NOT attempt to charge non-rechargeable batteries.
 - o DO NOT autoclave or spray batteries, battery contacts, charger, or AC power adapter with liquid of any kind. If corrosion appears on the contacts of the battery charger, call Ultradent Customer Service to order a replacement.
 - o DO NOT charge batteries around flammable materials.
 - o DO NOT keep charger in clinical operator.
- To avoid the risk of injury, DO NOT use batteries that are corroded (rust), dented, emit an odor or fluids, have a torn or missing wrapping, or are otherwise damaged. Call Ultradent Customer Service to order replacement batteries.
- To prevent the risk of thermal irritation or injury, avoid back-to-back curing cycles and do not expose oral soft tissues at close proximity for more than 10 seconds in any mode. If longer curing times are required, use multiple shorter curing cycles or use a dual-cure product to avoid heating soft tissue.
- Use caution when treating patients who suffer from adverse photobiological reactions or sensitivities, patients who are undergoing chemotherapy treatment, or patients being treated with photosensitizing medication.
- This unit may be susceptible to strong magnetic or static electric fields, which could disrupt the programming. If you suspect this has occurred, unplug the unit momentarily and then re-plug it into the outlet.

- DO NOT wipe down the curing light with caustic or abrasive cleaners, autoclave, or immerse in any kind of ultrasonic bath, disinfectant, cleaning solution, or liquid. Failure to follow included processing instructions may render curing light inoperable.
- To avoid damaging the equipment, DO NOT insert fingers, instruments, or other objects into the battery compartment of the curing light.
- To avoid damaging the equipment, DO NOT attempt to clean the gold contacts, or any part of the battery compartment. Call Ultradent Customer Service if there is a concern.
- To prevent the risk of cross-contamination, barrier sleeves are single patient use.
- To reduce the risk of corrosion, remove barrier sleeve after use.
- To reduce the risk of under-cured resins, do not use curing light if lens is damaged.

4. Stepwise Instructions

Preparation

1. Charge batteries before using the curing light (See section Battery Maintenance).
2. Place the curing light into a standard dental unit mounting bracket or accessory mounting bracket until ready for use.
3. Prior to each use, place a new barrier sleeve over the curing light.

Installing Hygienic Barrier Sleeves:

The hygienic barrier sleeve is custom fitted to the curing light and keeps the surface of the curing light clean. The barrier sleeve helps prevent cross contamination, helps keep dental composite material from adhering to the surface of the lens and curing light, and prevents discoloration and corrosion from cleaning solutions.

Note:

- Using the hygienic barrier sleeve will reduce the light output by 5-10%. Due to the high output power of the curing light, curing has been shown to be substantially equivalent.
- The curing light must be cleaned and sanitized with appropriate cleaning and/or sanitizing agents after each patient. See section titled Processing.

VALO Cordless Light Shield:

The VALO Cordless Light Shield is oval-shaped, can be rotated for maximum use, and can be used with a transparent barrier sleeve.

Use

1. Each power mode is used for the curing of dental materials with photo initiators. See Quick Mode Guide for recommended curing times.

NOTE: The curing light is programmed to cycle from the Standard Power to the High Power to the Xtra Power mode in sequence. For example, to change from the Standard Power mode to the Xtra Power mode, it is necessary to cycle into the High-Power mode and then to the Xtra Power mode.

2. The curing light stores the most recently used timing interval and mode, and it will default back to this whenever the modes are changed or if the batteries are removed.

Operation

CURING MODE: Standard Power mode

TIMING INTERVALS: 5, 10, 15, 20 seconds.

- The curing light defaults to this mode when it is INITIALLY powered on. The Mode/Status Light will be green and the four green Timing Lights are illuminated, indicating Standard Power mode.
- To change timing intervals, quickly press the Time/Mode Button.
- Press the Power Button to cure. To stop curing prior to completion of a timing interval, press the Power Button again.

CURING MODE: High Power mode

TIMING INTERVALS: 1, 2, 3, 4 seconds.

- From Standard Power mode, press and hold the Time/Mode Change Button for 2 seconds and release. The Mode/Status Light will be orange, and the four green Timing Lights will illuminate and flash, indicating High Power mode.
- To change timing intervals, quickly press the Time/Mode Button.
- Press the Power Button to cure. To stop curing prior to the completion of a timing interval, press the Power Button again.
- To return to Standard Power mode, press and hold the Time/Mode Change Button for 2 seconds and release, this will cycle to Xtra Power mode. Press and hold again for 2 seconds, and release. The Mode/Status Light will be green and the four green Timing Lights are illuminated, indicating Standard mode.

CURING MODE: Xtra Power mode

TIMING INTERVAL: 3 seconds only (Note: The Xtra Power mode has a 2 second safety delay at the end of each curing cycle to limit heating during consecutive curing. At the end of the delay, beeping indicates unit is ready for continued use).

- From Standard Power mode, press the Time/Mode Change Button for 2 seconds, release, press and hold again for 2 seconds, and release. The Mode/Status Light will be orange and flash, and three of the green Timing Lights will illuminate and flash, indicating Xtra Power mode.
- Press the Power Button to cure. To stop curing prior to the completion of a timing interval, press the Power Button again.
- To return to the Standard Power mode, press and hold the Time/Mode Button for 2 seconds and release. The Mode/Status Light will be green and the green Timing Lights are illuminated, indicating Standard Power mode.

Sleep Mode: The curing light will go into Sleep Mode after 60 seconds of inactivity, as indicated by a slow flashing of the mode/status light. Picking up or touching the unit will wake-up the curing light and automatically return it to the last setting used. To maximize battery life, leave the curing light undisturbed when not in use.

Cleanup

1. Discard used barrier sleeves in standard waste after each patient.
2. See Processing Section.

Mounting Bracket Instructions

1. Bracket should be mounted to a flat, oil-free surface.
2. Clean surface with rubbing alcohol.
3. Peel backing off the bracket's adhesive tape.
4. Position bracket so the curing light lifts upward when removed. Press firmly into place.

Quick Mode Guide

Mode	Standard Power				High Power				Xtra Power
Power Button									
Mode/Timing LEDs									
Time Buttons									
Time Options	5s	10s	15s	20s	1s	2s	3s	4s	3s Only
To Change Time	Press and release Time Button quickly to cycle through time options.								
To Change Modes	Press and hold Time Button for 2 seconds and release. VALO will cycle to next Mode.								
Legend	Solid LEDs				Blinking LEDs				

Quick Curing Guide

Recommended Curing Times for Optimal Results with VALO			
Mode	Standard Mode	High Power Mode	Xtra Power Mode
Per Layer	One 10 second cure	Two 4 second cures	One 3 second cure
Final Cure	Two 10 second cures	Three 4 second cures	Two 3 second cures
<p>Note: Exposure settings and times may need to be adjusted due to composite reactivity, shade, distance from the light lens to the composite, and depth of composite layer. It is up to the dental professional to know the requirement for the material they are using to determine the adequate time and settings.</p>			

Quick Warning Guide:

Power Level Warning	Temperature Warning	Service Warning	LED Warning
Replace batteries	Allow for cool down	Call Customer Service for Repair	Call Customer Service for Repair
<ul style="list-style-type: none"> • Low battery: slow flashing • Shut off: 3 beeps, flashing • Prohibits operation 	<ul style="list-style-type: none"> • 3 beeps • Slow flashing • Prohibits operation 	<ul style="list-style-type: none"> • No sound • Flashing, 2 seconds • Allows operation 	<ul style="list-style-type: none"> • Continuous 3 beeps • Single rapid flash • Prohibits operations

5. Maintenance Repair

User-Performed Repair

1. Routinely check the lens for cured dental resins. If necessary, use a non-diamond dental instrument to carefully remove any adhered resin.
2. Light meters differ greatly and are designed for specific light guide tips and lenses. Ultradent recommends routinely checking the output in Standard Power mode. NOTE: the true numeric output will be skewed due to the inaccuracy of common light meters and the custom LED pack in the curing light.

Manufacturer Repair

1. Repairs are only to be performed by authorized service personnel. Ultradent to provide service personnel with documentation to perform repairs.
2. When sending units in for repair, service, or calibrations, always remove the batteries from the curing light and charger. Wrap batteries, charger, adapter, and the curing light separately in the return box.
3. Ship batteries in accordance with local regulations.

Battery Maintenance

Charging and Charging Batteries

The curing light comes with 4 rechargeable lithium iron phosphate batteries.

How to charge batteries:

1. Plug charger into electrical outlet.
2. Insert batteries into charger with positive (+) end pointed towards the indicator lights on charger.
3. Green Lights indicate batteries are ready for use.
4. Batteries will take 1-3 hours to charge. Leave batteries in charger until ready for use.

NOTE: If the red light on the charger does not turn green when the batteries have been charging longer than three hours, the battery may be degraded (corrupted) and cannot be charged. Try a new battery or call Ultradent Customer Service to order a new set of rechargeable batteries.

How to change/insert batteries:

1. Remove back cap by twisting counterclockwise one quarter turn.
2. Remove batteries.
3. Insert fresh batteries positive (+) side first.
4. Reattach back cap by aligning and gently pushing while twisting clockwise. The cap will click when fully attached.
5. The unit is ready for use.

If necessary, Ultradent authorizes the following NON-RECHARGEABLE batteries for the curing light:

NOTE: DO NOT attempt to charge non-rechargeable batteries.

1. Tenergy Propel Photo Lithium
2. Titanium Innovations CR123A
3. Energizer® 123
4. Duracell® Ultra CR123A
5. SureFire® SF123A
6. Panasonic® CR123A

Low Batteries: The curing light signals the user that it is time to change the batteries when the Low Battery Indicator Light is flashing red. If the battery charge becomes too low, an audible 3 beep warning sound will occur and the curing light will not allow further operation until batteries are recharged or new batteries are inserted. (See Quick Warning Guide)

Charging Time and Battery Life: Fully charged battery life in the curing light is dependent on the Mode/Time interval, battery type, amount of use, and LED efficiency. In general, rechargeable batteries should last 1 - 2 weeks. Non-rechargeable batteries may last 2 - 3 times longer.

- Recommended Recharge Interval: When low battery indicator comes on, or approximately every 1 – 2 weeks, depending on use.
- Extra Batteries: The curing light comes with 4 rechargeable batteries. We also recommend keeping a spare set of non-rechargeable CR123A Batteries on hand in the event of a problem or loss of the rechargeable batteries.
- Battery Life Expectancy: Rechargeable Lithium Iron Phosphate batteries can be recharged approximately 1000 to 2000 times. Based on normal use and proper care, batteries should last up to five years, but should be replaced as necessary.

Warranty

Ultradent hereby warrants that this instrument shall, for a period of 5 years*, conform in all material respects to the specifications therefore as set forth in Ultradent's documentation accompanying the product and be free from any defects in materials or workmanship. This warranty applies solely to the original purchaser and is not transferable. All defective products are to be returned to Ultradent. There are no user service components of the VALO Cordless curing light system. Tampering with the VALO Cordless curing light will void its warranty. The VALO Cordless curing light warranty does not cover customer damage. For example; if a VALO Cordless is misused or dropped and the lens breaks, the customer would be responsible to pay for any necessary repairs.

*With sales receipt indicating the date of sale to the dentist.

6. Processing

After each use, moisten a gauze or soft cloth with an approved surface disinfectant and wipe the surface and lens.

ACCEPTABLE CLEANERS:

- Lysol Brand III Disinfectant Spray (Recommended)
- Isopropyl alcohol
- Ethyl alcohol-based cleaners
- Lysol™ Concentrate (alcohol-based only)

UNACCEPTABLE CLEANERS - DO NOT USE:

- Strong alkali detergent of any type, including hand soaps and dish soaps
- Bleach-based cleaners (e.g. Clorox™, Sterilox™)**
- Hydrogen Peroxide based cleaners
- Abrasive Cleaners (e.g. Comet Cleanser™)**
- Acetone or hydrocarbon-based cleaners
- MEK (Methyl Ethyl Ketone)
- Birex™**
- Gluteraldehyde
- Quaternary Ammonium Chloride salt-based cleaners
- Cavicide1™** solution or wipes
- Cavicide™** products (Non-Bleach)**

*Trademark of a company other than Ultradent

** If used, it may fade the color

BATTERY CHARGER:

If cleaning becomes necessary then unplug charger, dampen a cloth with isopropyl alcohol, and gently wipe the surface of the charger or batteries. Let the charger dry completely before resuming use.

CLEANING THE LIGHT SHIELD:

Cold disinfect the VALO Cordless Light Shield using any surface disinfectant. DO NOT autoclave.

7. Storage and Disposal

If storing the curing light for periods longer than 2 weeks, or packing it for travel, always remove the batteries. If batteries are left in the unit for long periods of time without recharging they may become nonfunctional or un-chargeable. Do not store batteries in temperatures over 60°C (140°F) or in direct sunlight.

Curing light Storage and Transport:

- Temperature: +10°C to +40°C (+50°F to +104°F)
- Relative Humidity: 10% to 95%
- Ambient Pressure: 500 hPa to 1060 hPa

When disposing of electronic waste, (i.e. devices, chargers, batteries and power supplies), follow local waste and recycling guidelines.


8. Technical Considerations

Accessories

Item	CE Information		
VALO Barrier Sleeves	  MDSS GmbH Schiffgraben 41 30175 Hanover Germany	Manufactured by: TIDI Products, LLC 570 Enterprise Drive Neenah, WI 54956 Made in USA	Distributed by: Ultradent Products Inc 505 West Ultradent Drive (10200 South) South Jordan, UT 84095 USA
VALO Cordless Light Shield			
VALO Charger			
VALO Batteries			

Technical information/Data

Effective Composite Curing Wavelength Bands:

Attribute	Information/Specification					
Lens	Diameter 9.75 mm					
Wavelength range	<ul style="list-style-type: none"> Utilizable wavelength range: 385 – 515nm Peak wavelengths: 395 – 415nm and 440 – 480nm 					
Light Intensity Table	Nominal Radiant Exitance Comparison Chart †				Radiant Exitance will vary based on instrument capability, measurement method and light placement. † Demetron radiometers and MARC spectrum analyzers should be used as reference only due to having smaller apertures than VALO curing lights. * Demetron radiometers should be used as reference only due to limitations in power and spectral response. ‡ Radiant Exitance conforms to ISO 10650 when measured with a Gigahertz spectrum analyzer.	
	Measurement Instrument	†* Demetron L.E.D. Radiometer	† MARC spectrum analyzer	‡ Gigahertz spectrum analyzer		
	Aperture of Meter	7 mm	3.9 mm	15 mm		15 mm
	Standard Power (±10%)	1000 mW/cm ²		900 mW/cm ²		670 mW
	High Power (±10%)	1400 mW/cm ²		1300 mW/cm ²		970 mW
	Xtra Power (±10%)		3200 mW/cm ² (+/-20%)	2100 mW/cm ²		1570 mW
VALO Cordless Curing Light	Ratings: IEC 60601-1 (Safety), IEC 60601-1-2 (EMC)		Weight: <ul style="list-style-type: none"> With batteries: 6 oz. (190 grams) Without batteries: 5 oz. (150 grams) Dimension: (8 x 1.28 x 1.06) inches, (203 x 32.5 x 27) mm			
Charger Power Supply	Output - 12VDC at 500mA Input - 100VAC to 240VAC Ultraviolet P/N 5930 VALO Charger Power Supply with international plug inserts		Ratings: IEC 60601-1 (Safety) Cord Length - 6 feet (1.8 meters) VALO Charger Power Supply is a Medical Grade Class II power supply and provides isolation from MAINS power			
VALO Charger	VALO 3.6VDC Lithium Iron Phosphate smart battery charger: <ul style="list-style-type: none"> Automatic shut off when fully charged Auto-detection of defective batteries Protections: Thermal, Overcharge, Short-circuit, reverse polarity <ul style="list-style-type: none"> Red LED – Charging Green LED – Empty or Fully Charged LED off – short circuit Charging time: 1 – 3 hours Rating: CE, WEEE					
VALO Batteries	Rechargeable: Safe chemistry Lithium Iron Phosphate (LiFePO4) RCR123A <ul style="list-style-type: none"> Working Voltage: 3.2VDC Ratings: UL, CE, RoHS, WEEE					
Operating Conditions	Temperature: +10°C to +32°C (+50°F to +90°F) Relative Humidity: 10% to 95% Ambient Pressure: 700 hPa to 1060 hPa					
Duty Cycle:	The curing light is designed for short-term operation. At maximum ambient temperature (32°C) 1 minute ON back-to-back cycling, 30 minutes OFF (cooling-off period).					

Trouble Shooting

If the solutions suggested below do not rectify the problem, please call Ultradent at 800.552.5512. Outside the United States, call your Ultradent distributor or dental dealer.	
Problem	Possible Solutions
Light will not turn on	<ol style="list-style-type: none"> 1) Press the Time/Mode or Power Button to wake from Power Save Mode. 2) Check the red Low Battery Indicator for battery charge status. 3) Check that fresh batteries are correctly inserted into the unit. 4) If red and yellow Warning LEDs are flashing this means the curing light has reached its internal temperature safety limit. Allow the curing light to cool down for 10 minutes or use a cool moist towel to cool the unit down quickly. 5) If red Warning LED flashes and beeps continuously, call Ultradent Customer Service for repair.
Light does not stay on for desired time	<ol style="list-style-type: none"> 1) Check Mode and Timing lights for correct time input. 2) Check the Low Battery Indicator for battery charge status. 3) Check that fresh batteries are properly inserted into the unit.
Light is not curing resins properly	<ol style="list-style-type: none"> 1) Check lens for residual cured resins/composites. 2) Using proper amber UV eye protection, verify the LED lights are working. 3) Check power level with light meter. If using a light meter, Ultradent recommends checking the curing light in Standard Power mode. <p>NOTE: The true numeric output will be skewed due to the inaccuracy of common light meters and the custom LED pack the curing light uses. Light meters differ greatly and are designed for specific light guide tips and lenses.</p> <ol style="list-style-type: none"> 4) Check expiration date on curing resin. 5) Ensure proper technique is being followed according to manufacturer recommendations.
Batteries will not charge	<ol style="list-style-type: none"> 1) Make sure batteries are inserted into the charger in the correct orientation and allow batteries to charge for 1-3 hours. 2) If red lights on the charger do not change to green, call Ultradent Customer Service to order replacement batteries and/or charger. 3) If neither green nor red lights on the charger are visible, call Ultradent Customer Service to order or replace charger and/or AC adapter.
Charger does not charge batteries	<ol style="list-style-type: none"> 1) Make sure charger is plugged in and AC adapter is plugged into a working power outlet. 2) If green or red lights on the charger are not visible, call Ultradent Customer Service for new charger and/or AC adapter.
Cannot change mode or time intervals	<ol style="list-style-type: none"> 1) Hold both Time/Mode and Power buttons down until a series of beeps indicates the curing light is unlocked.

9. Miscellaneous Information

Guidance and Manufacturer's Declaration for Electromagnetic Emissions		
The curing light is intended for use in the electromagnetic environment specified below. The customer or user should ensure that it is used in such an environment.		
Emissions Test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The curing light uses batteries and is not affected by EMI, RF, or surge suppression.
RF emission CISPR 11	Class B	The curing light uses electrical and electromagnetic energy only for their internal functions. Therefore, any RF emissions are very low and are not likely to cause interference in nearby electronic equipment.
Harmonic emissions IEC 61000-3-2	N/A	Harmonic emission and voltage fluctuation testing are not applicable to the curing light because it is battery powered.
Voltage fluctuations/flicker emissions IEC 61000-3-3	N/A	The curing light is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings for domestic use.

Guidance and Manufacturer's Declaration for Electromagnetic Immunity			
The curing light is intended for use in the electromagnetic environment specified below. The customer or user should ensure that it is used in such an environment.			
IMMUNITY test	IEC 60601 test level	Compliance level	Electromagnetic environment guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 8 kV contact ± 15 kV air	± 8 kV contact ± 15 kV air	Physical environment should be restricted to the following: 1. IP Code: IP20 2. Do not immerse in liquid. 3. Do not use around flammable gas. Unit is Non-APG and Non-AP. 4. Storage humidity range: 10% - 95% 5. Storage temperatures range: 10° C - 40° C
Electrical fast transient/ burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	± 2 kV for power supply lines Note 1: The curing light has no I/O ports	The curing light is battery powered and is not capable of connection to AC MAINS power.
Surge IEC 61000-4-5	± 1 kV line to line ± 2 kV line to earth	± 1 kV line to line ± 2 kV line to earth	Because the curing light is battery powered, it is not subject to electric transients, surges, voltage dips, shorts, interruptions, or variations on AC MAINS power.
Voltage, dips, shorts, interruptions and variations on the power supply input lines IEC 61000-4-11	<5% U (>95% dip in U for 0.5 cycle) 40% U (60% dip in U for 5 cycles) 70% U (30% dip in U for 25 cycles) <5% U (>95% dip in U for 5 s)	<5% U (>95% dip in U for 0.5 cycle) 40% U (60% dip in U for 5 cycles) 70% U (30% dip in U for 25 cycles) <5% U (>95% dip in U for 5 s) Note 2: Self recovers	The accessory battery charger may be subject to the above but it is separate from and not critical to the operation of the VALO Cordless curing light. If battery voltage of the curing light drops to 4VDC the unit will not allow operation. The VALO Cordless curing light will turn off. When new batteries are put in and the proper power levels are restored, the curing light will restart and return to the same state before power loss. The curing light will self-recover in the event of power loss.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical, residential, home health care, commercial, hospital, or military environment.
Note 1: The curing light is not equipped with any ports or any accessible I/O lines.			

Guidance and Manufacturer's Declaration for Electromagnetic Immunity for non-life support systems			
The curing light is intended for use in the electromagnetic environment specified below. The customer or user should ensure that it is used in such an environment.			
IMMUNITY test	IEC 60601 test level	Compliance level	Electromagnetic environment guidance
Conduction RF	3 Vrms	3 Vrms	Portable and mobile RF communications equipment should be used no closer to any part of the curing light, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance: $d = \left[\frac{3,5}{V_1} \right] \sqrt{P}$ $d = \left[\frac{3,5}{E_1} \right] \sqrt{P} \quad 80 \text{ MHz to } 800 \text{ MHz}$ $d = \left[\frac{7}{E_1} \right] \sqrt{P} \quad 800 \text{ MHz to } 2,5 \text{ GHz}$ P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol: 
IEC 61000-4-6	150 kHz to 80 MHz	150 kHz to 80 MHz	
Radiated RF	3 V/m	3 V/m	
IEC 61000-4-3	80 MHz to 2.5 GHz	80 MHz to 2.5 GHz	
NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies. NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people. a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the curing light is used exceeds the applicable RF compliance level above, the VALO Cordless curing light should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the curing light. b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.			

Guidance and Manufacturer's Declaration for recommended separation distances between portable and mobile RF communications equipment and the VALO Cordless curing light

The curing light is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The user of the curing light can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the curing light as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter (P in Watts)	Separation distance according to frequency of transmitter (meters)		
	150 kHz – 80 MHz $d = \left[\frac{3,5}{V_1} \right] \sqrt{P}$	80 MHz – 800 MHz $d = \left[\frac{3,5}{E_1} \right] \sqrt{P}$	800 MHz – 2.5 GHz $d = \left[\frac{7}{E_1} \right] \sqrt{P}$
0.01	0.12 meters	0.035 meters	0.07 meters
0.1	0.37 meters	0.11 meters	0.22 meters
1	1.7 meters	0.35 meters	0.7 meters
10	3.7 meters	1.11 meters	2.22 meters
100	11.7 meters	3.5 meters	7.0 meters

The curing light has been tested according to IEC 60601-1-2:2014 and passed under radiated field strengths of 10 V/m between 80-MHz to 2.5 GHz. The value of 3Vrms corresponds to V1 and the value 10 V/m corresponds to E1 in the formulas above.

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

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© 2021 Ultradent Products, Inc.
Ultradent Products, Inc.
Ultradent Products GmbH