## Peak™ Universal Bond and Peak™ SE

### Description:

**Peak Universal Bond** is a bottle or single syringe delivered bonding resin. It can be used with self-etch or with total etch. It is 7.5% filled with an ethyl alcohol solvent carrier and will cure with most high intensity lights including LEDs. Peak Universal Bond contains 0.2% chlorhexidine which may ensure long term bond strengths.

**Peak SE** Primer is a self-etching dentin bonding resin which is comprised of an etchant and primer in a syringe delivery system.

**Ultra-Etch™** is a syringe-delivered, 35% phosphoric acid etch and rinse gel.

Inclusion of the chlorhexidine additive has not shown to correlate with a reduction in secondary decay in patients. Clinical studies to evaluate the effects of chlorhexidine on oral bacterial or caries have not been performed.

### Indications for Use:

Use for virtually all etching and bonding needs in restorative dentistry for patients of all ages.


Peak Universal Bond bonds to the following materials:

- Dentin and Enamel
- Porcelain, zirconia
- Metal
- Composite

### Directions for Use:

**NOTE: Initial use of Peak SE requires activation of the two components separated in the JetMix™ syringe.**

1. Peak SE is stable for 120 days after activation.
2. Write the expiration date on the flat area of the plunger stem before activation.
3. DO NOT REMOVE LUER LOCK CAP PRIOR TO MIXING
4. Quickly and forcefully press "dated" plunger into middle barrel.
5. Verify that plunger is flush with middle barrel.

### Direct Bonding Technique and Porcelain/Zirconia Repair

#### 1. Prepare

- a. Dentin/Enamel: Isolate, clean and remove caries ensuring all non-mineral dentin is removed (Seek/Sable Seek is recommended for this). For abrasion/abfraction Class V preps, roughen with a diamond bur.
- b. Composite: Remove weakened areas of existing composite. Roughen surface with diamond bur.
- c. Metal: Microabrade surface.
- d. Zirconia: Clean surface and remove weakened zirconia with a diamond bur.
- e. Porcelain: Roughen and remove weakened porcelain with diamond bur. Clean fractured porcelain area for 5 seconds with Ultra-Etch. Rinse and dry.

#### 2. Etch

- a. Self-Etch Technique – dentin/enamel, composite and metal
  - i. Rinse preparation and leave damp.
  - ii. OPTIONAL – Apply Consepis™ (where available) with Black Mini™ Brush Tip to preparation with a scrubbing motion and lightly air

dry leaving surface damp.

- iii. Apply Peak SE with Black Mini Brush Tip to the tooth preparation for 20 seconds using continuous agitation on dentin. Do not scrub the enamel.
  - iv. Thin/dry for 3 seconds using the air/water syringe or by placing the high volume suction directly over the preparation.
  - v. Proceed to "Bond" step below.
- b. Total Etch Technique – dentin/enamel, composite and metal
- i. Verify the color and consistency of the product flow before applying by expressing a small amount onto a mixing pad or gauze. Discard any clear liquid that is expressed.
  - ii. Apply Ultra-Etch with the Blue Micro™ Tip to all surfaces of the tooth preparation for 20 seconds.
  - iii. Rinse thoroughly for 5 seconds.
  - iv. Lightly dry using the air/water syringe or by placing the high volume suction directly over the preparation. Leave the surface slightly damp.
  - v. OPTIONAL – Apply Consepis (where available) with Black Mini Brush Tip to preparation with a scrubbing motion and lightly air dry leaving surface damp.
  - vi. Proceed to "Bond" step below.
- c. Porcelain Etch – CAUTION – use proper isolation such as an acid neutralizing barrier (EtchArrest™) when using hydrofluoric acid (HF) (Ultradent™ Porcelain Etch). Do not allow hydrofluoric acid to touch gingiva or dentin.
- i. Apply hydrofluoric acid to fractured area for 90 seconds.
  - ii. Suction area from surface THEN thoroughly rinse and dry.
  - iii. Apply Ultra-Etch for 5 seconds to remove porcelain salts and debris formed by HF etching. Rinse and dry.
  - iv. Apply puddle coat of silane for 60 seconds.
  - v. Dry thoroughly. DO NOT RINSE.
  - vi. Proceed to "Bond" step below.
- d. Zirconia – DO NOT ETCH OR USE SILANE
- i. Air-abrade fractured surface of prosthesis, rinse and dry.
  - ii. Apply zirconia primer to fractured area following manufacturer's instructions.
  - iii. Proceed to "Bond" step below.
3. **Bond**
- a. Bottle Delivery
    - i. Dispense 1-3 drops of Peak™ Universal Bond into the "B" side of the disposable dappen dish.
    - ii. Using the Micro Applicator Brush, apply a puddle coat of Peak Universal Bond to the preparation and gently agitate for 10 seconds.
    - iii. Thin/dry for 10 seconds using ¼ to ½ air pressure. Preparation should appear shiny.
    - iv. Light cure for 10 seconds (20 seconds for lights with output <60mW/cm²).
    - v. Restore with flowable or packable composite as per manufacturer's instructions.
  - b. Syringe Delivery
    - i. Apply a puddle coat of Peak Universal Bond with Inspiral™ Brush Tip – gently agitate for 10 seconds.
    - ii. Thin/dry 10 seconds using ¼ to ½ air pressure. Preparation should appear shiny.
    - iii. Light cure for 10 seconds (20 seconds for lights with output <60mW/cm²).
    - iv. Restore with flowable or packable composite as per manufacturer's instructions.

## Indirect Bonding Technique - Indirect Crowns, In-lays, Onlays,

### Zirconia

#### 1. Prepare

- Remove temporary, clean preparation, rinse, and dry (Consepsis Scrub with a STARbrush™ to clean, disinfect, and remove residual cement).
- Verify prosthetic fit.
- Prepare inside surface of prosthesis.
  - Metal based:
    - Microabrade (sandblast) inside surface of prosthesis.
    - Rinse and dry surface.
  - Ceramic/porcelain:
    - Apply hydrofluoric acid (HF) (Ultradent Porcelain Etch) to inside surface of prosthesis for 90 seconds, rinse and dry.
    - Apply Ultra-Etch for 5 seconds to remove porcelain salts and debris formed by HF etching. Rinse and dry.
    - Apply silane (Ultradent™ Silane) to inside surface of prosthesis for 1 minute, dry and set prosthesis aside. Do not rinse.
  - Zirconia: DO NOT ETCH OR USE SILANE
    - Air abrade internal surface of prosthesis, rinse and dry.
    - Apply a zirconia primer as per manufacturer instructions.

#### 2. Etch - Prepare tooth surface

- Self-Etch Technique – dentin/enamel, composite and metal
  - Rinse preparation and leave damp.
  - OPTIONAL – Apply Consepsis (where available) with Black Mini Brush Tip to preparation with a scrubbing motion and lightly air dry leaving surface damp.
  - Apply Peak™ SE with Black Mini™ Brush Tip - 20 seconds to tooth preparation using continuous agitation on dentin. Do not scrub the enamel.
  - Thin/dry for 3 seconds using the air/water syringe or by placing the high volume suction directly over the preparation.
  - Proceed to “Bond” step below
- Total Etch Technique – dentin/enamel, composite and metal
  - Apply Ultra-Etch with the Blue Micro Tip to all surfaces of the tooth preparation for 20 seconds.
  - Rinse thoroughly for 5 seconds.
  - Lightly dry using the air/water syringe or by placing the high volume suction directly over the preparation. Leave the surface slightly damp.
  - OPTIONAL – Apply Consepsis (where available) with Black Mini Brush Tip to preparation with a scrubbing motion and lightly air dry leaving surface damp.
  - Proceed to “Bond” step below.
- Porcelain Etch – CAUTION – use proper isolation such as an acid neutralizing barrier (EtchArrest) when using hydrofluoric acid (HF) (Ultradent Porcelain Etch). Do not allow hydrofluoric acid to touch gingiva or dentin.
  - Apply hydrofluoric acid to fractured area for 90 seconds.
  - Suction acid from surface THEN thoroughly rinse and dry.
  - Apply Ultra-Etch for 5 seconds to remove porcelain salts and debris formed by HF etching. Rinse and dry.
  - Apply puddle coat of silane for 60 seconds.
  - Dry thoroughly. DO NOT RINSE.
  - Proceed to “Bond” step below.
- Zirconia – DO NOT ETCH OR USE SILANE
  - Air-abrade internal surface of prosthesis, rinse and dry.
  - Apply zirconia primer to fractured area following manufacturer’s instructions.
  - Proceed to “Bond” step below.

#### 3. Bond

- Bottle Delivery

- Dispense 1-3 drops of Peak Universal Bond into the “B” side of the disposable dappen dish.
  - Using the Micro Applicator Brush, apply a puddle coat of Peak Universal Bond to the preparation and gently agitate for 10 seconds.
  - Thin/dry for 10 seconds using full air pressure. Preparation should appear shiny.
  - Light cure for 10 seconds (20 seconds for lights with output <600mW/cm²).
- Syringe Delivery
    - Apply a puddle coat of Peak Universal Bond with Inspirial Brush Tip – gently agitate for 10 seconds.
    - Thin/dry 10 seconds using full air pressure. Preparation should appear shiny.
    - Light cure for 10 seconds (20 seconds for lights with output <600 mW/cm²).

#### 4. Cement

- Apply light cure cement for translucent veneers or dual cure resin cement to other prostheses and follow manufacturer’s instructions for use.

**NOTE:** Immediate Dentin Sealing requires Direct Technique steps 2A and 3 above. Apply glycerin over Peak Universal Bond and light cure 10 seconds. Rinse Glycerin. Make final impression. Re-apply glycerin to preparation as separating medium prior to cementing the temporary in place.

## Repair of Composite / Metal / Zirconia / Porcelain

- Composite:** Remove weakened areas of existing composite. Roughen surface with diamond bur.
  - Total-Etch Technique
    - Apply Ultra-Etch with the Blue Micro Tip to all surfaces of tooth preparation for 20 seconds.
    - Rinse thoroughly for 5 seconds and leave damp.
  - Metal:** Microabrade surface
    - Total-Etch Technique
      - Apply Ultra-Etch with the Blue Micro Tip to all surfaces of tooth preparation for 20 seconds.
      - Rinse thoroughly for 5 seconds and leave damp.
  - Zirconia:** Clean surface; remove weakened zirconia with diamond bur. DO NOT ETCH OR USE SILANE.
    - Air-abrade internal surface of prosthesis.
    - Apply zirconia primer to fractured area following manufacturer’s instructions.
  - Porcelain:** Roughen and remove weakened porcelain with diamond bur. Clean fractured porcelain area for 15 seconds with Ultra-Etch. Rinse and dry.

CAUTION – use proper isolation such as an acid neutralizing barrier (EtchArrest) when using hydrofluoric acid (HF) (Ultradent™ Porcelain Etch). Do not allow hydrofluoric acid to touch gingiva or dentin.

    - Apply hydrofluoric acid to fractured area for 90 seconds.
    - Suction acid from surface then thoroughly rinse and dry surface.
    - Apply Ultra-Etch for 5 seconds to remove porcelain salts and debris formed by HF etching. Rinse and dry.
    - Apply puddle coat of silane for 1 minute.
    - Dry thoroughly. DO NOT RINSE.
- Bond**
  - Bottle Delivery
    - Dispense 1-3 drops of Peak Universal Bond into the “B” side of the disposable dappen dish.
    - Using the Micro Applicator Brush, apply a puddle coat of Peak Universal Bond and gently agitate for 10 seconds.
    - Thin/dry for 10 seconds using 1/4 to 1/2 air pressure. Preparation

- should appear shiny.
- iv. Light cure for 10 seconds (20 seconds for lights with output <60mW/cm<sup>2</sup>).
- b. Syringe Delivery
  - i. Apply a puddle coat of Peak Universal Bond with Inspiral Brush Tip - gently agitate for 10 seconds.
  - ii. Thin/dry for 10 seconds using 1/4 to 1/2 air pressure.
  - iii. Preparation should appear shiny.
  - iv. Light cure for 10 seconds (20 seconds for lights with output <60mW/cm<sup>2</sup>).
- 6. Restore with flowable and/or packable composite as per manufacturer's instructions.

## Post and Core

### 1. Prepare

- a. Drill post hole using a rubber stop and appropriate size UniCore™ Drill.
- b. Thoroughly rinse the post space for 10 seconds using the TriAway Adaptor and a 22ga Endo-Eze™ Tip. Dry with TriAway Adaptor and 22ga Endo-Eze Tip from the bottom of the preparation up.
- c. Verify post (UniCore Post) fit and cut to desired length using a high speed diamond disc. NOTE: Carbide and operative burs may fray the fibers.

### 2. Etch

- a. Self-Etch Technique
  - i. Rinse canal space and leave damp.
  - ii. Attach a blue 22ga Deliver-Eze™ Brush Tip to Peak SE syringe.
  - iii. Apply to the post space and coronal preparation for 20 seconds using an agitating action.
  - iv. Blow out excess with TriAway Adaptor and 22ga Endo-Eze Tip from the bottom of the space, upward. Thin/dry for 3 seconds using the air/water syringe or by placing the high volume suction directly over the preparation.
  - v. Proceed to "Bond" step below.
- b. Total-Etch Technique:
  - i. Attach a blue 22ga Endo-Eze Tip to the Ultra-Etch syringe.
  - ii. Start apically and fill post space coronally. Apply to coronal preparation. Etch for 20 seconds.
  - iii. Suction off excess etchant.
  - iv. Using the TriAway Adaptor and 22ga Endo-Eze Tip, rinse thoroughly and lightly air dry, leaving the post space slightly damp.
  - v. Proceed to "Bond" step below.

### 3. Bond

- a. Bottle Delivery
  - i. Dispense 1-3 drops of Peak Universal Bond into the "B" side of the disposable dappen dish.
  - ii. Using the Nano-Brush™ or the NavTip FX, apply Peak Universal Bond to the post space and gently agitate the full length of the canal and tooth preparation for 10 seconds.
  - iii. Thoroughly remove excess using the TriAway Adaptor and the 22ga Endo-Eze Tip.
  - iv. Air dry post space and preparation for 10 seconds using full air pressure to thin/dry with TriAway Adaptor and 22ga Endo-Eze Tip to evaporate any remaining solvent from the bottom up.
  - v. Light cure down the post space for 20 seconds (40 seconds for lights with output <60mW/cm<sup>2</sup>).
  - vi. Use paper point or previously fitted post to verify post space is not blocked.
- b. Syringe Delivery
  - i. Attach a 30ga NavTip FX to the Peak Universal Bond syringe.
  - ii. Place 1-2 drops into post space and onto preparation, and then agitate full length of the canal, including the entire tooth preparation, for 10 seconds.

- iii. Thoroughly remove excess Peak Universal Bond using the TriAway Adaptor and the 22ga Endo-Eze Tip.
- iv. Air dry post space and preparation for 10 seconds using full air pressure to thin/dry with TriAway Adaptor and 22ga Endo-Eze Tip to evaporate any remaining solvent from the bottom up.
- v. Light cure down the post space for 20 seconds (40 seconds for lights with output of <60mW/cm<sup>2</sup>). Use paper point or previously fitted post to verify post space is not blocked.

## Post Cementation and Core Build-up

1. Remove cap from the PermaFlo™ DC syringe.
  2. Place a mixing tip onto the dual barrel syringe lining up the two ports.
  3. Twist clockwise until locked in place.
  4. Attach the flexible Intraoral Tip™ to the mixing tip for delivery directly into the canal from the dual barrel syringe.
  5. Express a small amount before placing in the post space.
  6. Insert the tip into the post space to full depth.
  7. Using even pressure deliver PermaFlo DC starting apically and moving coronally.
  8. Insert the (UniCore™ Post) slowly, displacing the excess cement.
  9. Place curing light directly above the post and light cure 20 seconds; 40 seconds for lights with output <60mW/cm<sup>2</sup>.
  10. Continue building layers of PermaFlo™ DC around post as needed. NOTE: Light cure for 5 seconds between layers to gel and form core. Do a final cure on core for 40 seconds.
- Do not remove mixing tip until next use.** Disinfect syringe and mixing tip.

## Precautions and warnings:

1. Carefully read and understand all instructions before using Peak Universal Bond.
2. For professional use only.
3. Re-cap immediately following use to avoid polymerization.
4. Always verify material flow of syringe materials prior to applying intraorally. If resistance is met, replace tip and re-check.
5. Keep caps on the bottle until use.
6. Clinician and patient should wear UV protective eyewear when curing resin materials.
7. If not used daily, refrigerate product.
8. Bring refrigerated products to room temperature before using.
9. To optimize bond strengths use oil-free and moisture-free air.
10. Redirect overhead light to prevent premature polymerization of all resin based materials.
11. Resins can be sensitizing. Avoid repeated contact of uncured dental resin with skin. Do not use on patients with a known sensitivity to acrylates or other resins. If allergic reaction, dermatitis or rash develops, consult a physician.
12. Dispose of tips, empty syringes and bottles properly.
13. Keep out of reach of children.
14. Do not use after expiration date noted on containers.
15. Re-cap, disinfect and wipe syringe with an intermediate level disinfectant between uses. If disposable syringe cover is used, remove tip, re-cap, and discard syringe cover.
16. Tips are disposable. To avoid cross-contamination, do not re-use tips.
17. Keep products out of direct heat/sunlight.
18. Isolate strong chemicals to area of treatment.

**NOTE: For SDS and additional information about using Peak Universal Bond, Peak SE, and other related products, please go to [www.ultradent.com](http://www.ultradent.com)**

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



### Peak™ Universal Bond


For product SDS please see our website: [www.ultradent.com](http://www.ultradent.com)

### Peak™SE

For product SDS please see our website: [www.ultradent.com](http://www.ultradent.com)

### Ultra-Etch™

For product SDS please see our website: [www.ultradent.com](http://www.ultradent.com)

#### HAZARD RATING

4 = Severe  
3 = Serious  
2 = Moderate  
1 = Slight  
0 = Minimal

U.S. federal law restricts the sale of this device by or on the order of a dentist. For immediate reorder and/or complete descriptions of Ultradent's product line, refer to Ultradent's catalog or call Toll Free **1-800-552-5512**.

Outside U.S. call (801) 572-4200  
or visit [www.ultradent.com](http://www.ultradent.com).

<sup>1</sup> Evaluation of experimental Peak containing chlorhexidine (CHX) blended at 0.2% within it's formulation. re: Professor Lorenzo Breschi from the Department of Biomedicine, University of Trieste, Trieste Italy.

**EC REP**

Ultradent Products GmbH  
Am Westthover Berg 30  
51149 Cologne  
Germany

© Copyright 2017 Ultradent Products, Inc. All Rights Reserved.

Manufactured by Ultradent Products, Inc.  
505 West Ultradent Drive (10200 South),  
South Jordan, Utah 84095, USA  
Made in USA

CE 0044

**ULTRADENT**  
PRODUCTS, INC.

45414.6 092419