# VALO™

ORTHO·CORDLESS

**ULTRADENT**  
PRODUCTS, INC.  
*Improving Oral Health Globally*


## 1. Product Description

With its broadband spectrum, the VALO Ortho Cordless curing light is designed to polymerize all light cured products in the wavelength range of 385—515nm per ISO 10650. The VALO Ortho Cordless curing light uses the Ultradent VALO rechargeable batteries and battery charger. The curing light is designed to rest in a standard dental unit bracket or can be custom mounted using the bracket included with the kit.

VALO Ortho Cordless Product Components:

- 1 – VALO Ortho Cordless curing light
- 4 – Ultradent VALO rechargeable batteries
- 1 – Ultradent VALO battery charger with medical grade 12VDC AC power adapter
- 1 – VALO™ Barrier Sleeve sample pack
- 1 – VALO™ Black Light lens
- 1 – VALO™ Cordless Light Shield
- 1 – Curing light surface mounting bracket with double stick adhesive tape

### Overview of Controls:

For all products described, carefully read and understand all instructions and SDS information prior to use.

## 2. Indications for Use/Intended Purpose

The source of illumination for curing photo-activated dental restorative materials and adhesives.

## 3. Warnings and Precautions

### Risk Group 2

CAUTION UV emitted from this product. Eye or skin irritation may result from exposure. Use appropriate shielding.

CAUTION Possibly hazardous optical radiation emitted from this product. Do not stare at operating lamp. May be harmful to the eyes.

- DO NOT look directly into the light output. Patient, clinician, and assistants should always wear amber colored UV eye protection when curing light is in use.
- To prevent the risk of electric shock, no modification of this equipment is allowed. Use only the included Ultradent VALO light power supply and plug adapters. If these components are damaged, do not use and call Ultradent Customer Service to order a replacement.
- Portable RF communications equipment may degrade performance if used closer than 30 cm (12 in.).
- Use only authorized accessories, cables, and power supplies to prevent improper operation, increased electromagnetic emissions or decreased electromagnetic immunity (refer to Electromagnetic Emissions section).
- To avoid the risk of electrical fire associated with handling of batteries:
  - o DO NOT mix rechargeable batteries with non-rechargeable batteries or other battery types.
  - o DO NOT attempt to charge non-rechargeable batteries. Charging non-rechargeable batteries can result in injury or property damage.
  - o DO NOT autoclave or spray batteries, battery contacts, charger, or AC power adapter with liquid of any kind. If corrosion appears on the contacts of the battery charger, call Ultradent Customer Service to order a replacement.
  - o DO NOT charge batteries around flammable materials.
  - o DO NOT keep charger in clinical operatory.
- To avoid the risk of injury, DO NOT use batteries that are corroded (rust), dented, emit an odor or fluids, have a torn or missing wrapping, or are otherwise damaged. Call Ultradent Customer Service to order replacement batteries.
- Only use recommended batteries. Alternative batteries can cause malfunction.

- To prevent the risk of thermal irritation or injury, avoid back-to-back curing cycles and do not expose oral soft tissues at close proximity for more than 10 seconds in any mode. If longer curing times are required, use multiple curing cycles with periods of rest between cycles, or use a dual-cure product to avoid heating soft tissue.
- Use caution when treating patients who suffer from adverse photobiological reactions or sensitivities, patients who are undergoing chemotherapy treatment, or patients being treated with photosensitizing medication.
- This unit may be susceptible to strong magnetic or static electric fields, which could disrupt the programming. If you suspect this has occurred, unplug the unit momentarily and then re-plug it into the outlet.
- DO NOT wipe down the curing light with caustic or abrasive cleaners, autoclave, or immerse in any kind of ultrasonic bath, disinfectant, cleaning solution, or liquid. Failure to follow included processing instructions may render curing light inoperable.
- To avoid damaging the equipment, DO NOT insert fingers, instruments, or other objects into the battery compartment of the curing light.
- To avoid damaging the equipment, DO NOT attempt to clean the gold contacts, or any part of the battery compartment. Call Ultradent Customer Service if there is a concern.
- To prevent the risk of cross-contamination, barrier sleeves are single-use.
- To reduce the risk of corrosion, remove barrier sleeve after use.
- To reduce the risk of under-cured resins, do not use curing light if lens is damaged.
- This unit is programmed for orthodontic use. If using the VALO Ortho Cordless curing light for any other dental procedures, DO NOT exceed 10 seconds in the standard power mode. For any other dental procedure: avoid overheating the soft tissue by using two 10 second cures instead of a full 20 second cure.

#### 4. Stepwise Instructions

##### Preparation

1. Charge batteries before using the curing light (See section Battery Maintenance).
2. Prior to each use, place a new barrier sleeve over the curing light, minimizing wrinkles over the lens for best results.

##### Note:

- To help prevent cross-contamination and help keep dental composite material from adhering to the surface of the lens and wand body, an Ultradent approved barrier sleeve must be used over the VALO curing light with each use. Barrier sleeves are intended for single-patient use.

##### VALO Cordless Light Shield:

The VALO Cordless Light Shield is oval-shaped, can be rotated for maximum use, and can be used with a transparent barrier sleeve.

##### Use

1. Each power mode is used for the curing of dental materials with photo initiators. See Quick Mode Guide for recommended curing times.
- NOTE: The curing light is programmed to cycle from the Standard Power to the Xtra Power mode to the Xtra Power Quadrant mode in sequence. For example, to change from the Standard Power mode to the Xtra Power Quadrant mode, it is necessary to cycle into the Xtra Power mode and then to the Xtra Power Quadrant mode.
2. The curing light stores the most recently used timing interval and mode, and it will default back to this whenever the modes are changed or if the batteries are removed.

##### Operation

###### CURING MODE: Standard Power mode

TIMING INTERVALS: 5, 10, 15, and 20 seconds.

- The curing light defaults to this mode when it is INITIALLY powered on. The Mode/Status Light will be green and the four green Timing Lights are illuminated, indicating Standard Power mode.
- To change timing intervals, quickly press the Time/Mode Button.
- Press the Power Button to cure. To stop curing prior to completion of a timing interval, press the Power Button again.

###### CURING MODE: Xtra Power mode

TIMING INTERVAL: 1, 2, and 3 seconds

- From Standard Power mode, press the Time/Mode Change Button for 2 seconds and release. The Mode/Status Light will be orange and flash, and three of the green Timing Lights will illuminate and flash, indicating Xtra Power mode.
- To change the timing intervals, quickly press the Time/Mode change button.
- Press the Power Button to cure. To stop curing prior to the completion of a timing interval, press the Power Button again.
- To return to the Standard Power mode, press and hold the Time/Mode Button for 2 seconds, release, hold for 2 seconds and release. The Mode/Status Light will be green and the green Timing Lights are illuminated, indicating Standard Power mode..

### CURING MODE: Xtra Power Quadrant mode

The Xtra Power Quadrant modes is designed to allow the clinician to continuously cure one quadrant, or five teeth without needing to remove the VALO Ortho Cordless hand piece from the oral cavity. Each of the five teeth will be exposed for 3 seconds followed by a 2 second rest.

**TIMING INTERVAL:** 5 automatically timed bursts of 3 seconds each (Note: there is a 2 second safety delay at the end of each curing cycle).

- From Standard Power mode, press the Time/Mode Change Button for 2 seconds, release, press and hold again for 2 seconds, and release. The Mode/Status Light will be orange and flash, and three of the green Timing Lights will illuminate and flash, indicating Xtra Power Quadrant mode.
- Press the Power Button to cure. To stop curing prior to the completion of a timing interval, press the Power Button again.
- To return to the Standard Power mode, press and hold the Time/Mode Button for 2 seconds and release. The Mode/Status Light will be green and the green Timing Lights are illuminated, indicating Standard Power mode.

**Sleep Mode:** The curing light will go into Sleep Mode after 60 seconds of inactivity, as indicated by a slow flashing of the mode/status light. Picking up or touching the unit will wake-up the curing light and automatically return it to the last setting used. To maximize battery life, leave the curing light undisturbed when not in use.

### Cleanup

- Discard used barrier sleeves in standard waste after each patient.
- See Processing Section.


### Mounting Bracket Instructions

- Bracket should be mounted to a flat, oil-free surface.
- Clean surface with rubbing alcohol.
- Peel backing off the bracket's adhesive tape.
- Position bracket so the curing light lifts upward when removed. Press firmly into place.

### Quick Mode Guide

SUGGESTED CURING TIMES	STANDARD	XTRA POWER	XTRA POWER QUADRANT
<b>POWER (mW/cm<sup>2</sup>)</b>	1200* 	3200 	3200 
<b>POWER BUTTON</b>	   	  	
<b>TIMING INTERVAL LIGHTS</b>			
<b>MODE/STATUS LIGHT</b>			
<b>TIME/MODE CHANGE BUTTON</b>	   	  	
	5 sec. 10 sec. 15 sec. 20 sec.	1 sec. 2 sec. 3 sec.	3 sec.
<b>METAL &amp; CERAMIC BRACKETS (CURE TIME PER BRACKET)</b>	1 X 10 SECONDS	2 X 3 SECONDS	2 X 3 SECONDS
			5 consecutive burst at 3 seconds each
<b>RESTORATIVE CURING (2MM INCREMENTS)</b>			
<b>PER LAYER</b>	1 X 10 SECONDS	1 X 3 SECONDS	—
<b>FINAL CURE</b>	1 X 20 SECONDS	2 X 3 SECONDS	—

### Quick Warning Guide:

Power Level Warning	Temperature Warning	Service Warning	LED Warning
Replace batteries	Allow for cool down	Call Customer Service for Repair	Call Customer Service for Repair
<ul style="list-style-type: none"> <li>Low battery; slow flashing</li> <li>Shut off; 3 beeps, flashing</li> <li>Prohibits operation</li> </ul>	<ul style="list-style-type: none"> <li>3 beeps</li> <li>Slow flashing</li> <li>Prohibits operation</li> </ul>	<ul style="list-style-type: none"> <li>No sound</li> <li>Flashing, 2 seconds</li> <li>Allows operation</li> </ul>	<ul style="list-style-type: none"> <li>Continuous 3 beeps</li> <li>Single rapid flash</li> <li>Prohibits operations</li> </ul>
  	  	  	  

## 5. Maintenance

### General Cleaning of Curing Light

After each use, moisten a gauze or soft cloth with an approved surface disinfectant and wipe the surface and lens. Unauthorized cleaners may cause damage to the curing light.

### ACCEPTABLE CLEANERS:

- 70% Isopropyl alcohol
- 70% Ethanol

### General Cleaning of the Light Shield:

Clean the VALO Cordless light shield using any surface disinfectant. DO NOT autoclave.

### Repair

#### User-Performed Maintenance

1. Use a barrier sleeve to prevent dental composite from adhering to the surface of the lens. If necessary, use a plastic or stainless steel dental instrument to carefully remove any adhered composite. Do not use tools that will damage the lens.
2. Light meters differ greatly and are designed for specific light guide tips and lenses. Ultradent recommends routinely checking the output in Standard Power mode. NOTE: the true numeric output will be skewed due to the inaccuracy of common light meters and the custom LED pack in the curing light.

### Battery Maintenance

#### Charging and Changing Batteries:

The curing light comes with 4 rechargeable lithium iron phosphate batteries. Ultradent recommends using VALO rechargeable batteries only.

#### How to charge batteries:

1. Place charger and batteries in an area where they will not come into contact with liquids or disinfectant spray.
2. Insert batteries into charger with positive (+) end pointed towards the indicator lights on charger.
3. Green Lights indicate batteries are ready for use.
4. Batteries will take 1–3 hours to charge. Leave batteries in charger until ready for use.

NOTE: If the red light on the charger does not turn green when the batteries have been charging longer than three hours, the battery may be degraded (corrupted) and cannot be charged. Try a new battery or call Ultradent Customer Service to order a new set of rechargeable batteries.

#### How to change/insert batteries:

1. Remove back cap by twisting counterclockwise one quarter turn.
2. Remove batteries.
3. Insert fresh batteries positive (+) side first.
4. Reattach back cap by aligning and gently pushing while twisting clockwise. The cap will click when fully attached.
5. The unit is ready for use.

Low Batteries: The curing light signals the user that it is time to change the batteries when the Low Battery Indicator Light is flashing red. If the battery charge becomes too low, an audible 3-beep warning sound will occur and the curing light will not allow further operation until batteries are recharged or new batteries are inserted. (See Quick Warning Guide)

Charging Time and Battery Life: Fully charged battery life in the curing light is dependent on the Mode/Time interval, battery type, amount of use, and LED efficiency. In general, rechargeable batteries should last about 400 cure cycles.

- Recommended Recharge Interval: When low battery indicator comes on, or approximately every 1–2 weeks, depending on use.
- Extra Batteries: The curing light comes with 4 rechargeable batteries.
- Battery Life Expectancy: Rechargeable Lithium Iron Phosphate batteries can be recharged approximately 1000 to 2000 times. Based on normal use and proper care, batteries should last up to five years, but should be replaced as necessary.

### Manufacturer Repair

1. Repairs are only to be performed by authorized service personnel. Ultradent is to provide service personnel with documentation to perform repairs.
2. When sending units in for repair, service, or calibrations always remove the batteries from the curing light and charger. Wrap batteries, charger, adapter, and the curing light separately in the return box.
3. Ship batteries in accordance with local regulations.

### Warranty

Ultradent Products, Inc. ("Ultradent") warrants that this product shall, for a period of 5 years from the date of purchase, when operated according to the operating instructions included with the product, (i) conform in all material respects to the specifications set forth in Ultradent's documentation accompanying the product; and (ii) be free from defects in material and workmanship.

This limited warranty is nontransferable and applies solely to the original purchaser and does not extend to subsequent owners of the product. This limited warranty does not cover any other accessory components such as, but not limited to, batteries, chargers, adapters, or adaptive lenses. This limited warranty is void if the product fails or is damaged due to negligence, abuse, misuse, accident, modification, tampering, alteration, or failure to follow the applicable instructions for use. For purposes of example only, a product that is dropped and damaged

is not covered under this warranty. To qualify under this limited warranty, proof of purchase (e.g., sales receipt or similar documentation) must be submitted to Ultradent along with the defective product.

A defective product meeting the warranty conditions set forth herein will, at Ultradent's sole discretion, either be repaired or replaced. In no event shall Ultradent's liability for the product exceed the purchase price paid by the purchaser. Under no circumstances shall Ultradent be liable for any indirect, incidental, foreseen, unforeseen, special, or consequential damages arising out of or in connection with the use of this product.

## 6. Processing

This product is intended to be used with a barrier sleeve. See above for maintenance and general cleaning.

## 7. Storage and Disposal

If storing the curing light for periods longer than 2 weeks, or packing it for travel, always remove the batteries. If batteries are left in the unit for long periods of time without recharging they may become nonfunctional or un-chargeable. Do not store batteries in temperatures over 40°C (104°F) or in direct sunlight.

Curing light Storage and Transport:

- Temperature: +10°C to +40°C (+50°F to +104°F)
- Relative Humidity: 10% to 95%
- Ambient Pressure: 500 hPa to 1060 hPa


When disposing of electronic waste, (i.e. devices, chargers, batteries and power supplies), follow local waste and recycling guidelines.

## 8. Technical Considerations

### Accessories

Item	CE Information
VALO Cordless Light Shield	
VALO Charger	
VALO Batteries	

### Technical information/Data

Attribute	Information/Specification		
Lens	Diameter 9.75 mm		
Wavelength range	<ul style="list-style-type: none"> <li>Utilizable wavelength range: 385 – 515nm</li> <li>Peak wavelengths: 395 – 415nm and 440 – 480nm</li> </ul>		
Light Intensity Table	Nominal Radiant Exitance Comparison Chart		
	Measurement Instrument	‡ Gigahertz-Optik spectrum analyzer	
		Exitance	Total Power
	Aperture of Meter	15 mm	15 mm
	Standard Power (±10%)	1100 mW/cm <sup>2</sup>	820 mW
	Xtra Power (±10%)	2100 mW/cm <sup>2</sup>	1570 mW
Xtra Power Quadrant (±10%)	2100 mW/cm <sup>2</sup>	1570 mW	
VALO Ortho Cordless Curing Light	Ratings: IEC 60601-1 (Safety), IEC 60601-1-2 (EMC)	Weight: <ul style="list-style-type: none"> <li>With batteries: 6 oz. (190 g)</li> <li>Without batteries: 5 oz. (150 g)</li> </ul> Dimension: (8 x 1.28 x 1.06) inches. (203 x 32.5 x 27) mm	
Charger Power Supply	Output - 12VDC at 500mA Input - 100VAC to 240VAC Ultradent P/N 5930 VALO Charger Power Supply with international plug inserts	Ratings: IEC 60601-1 (Safety) Cord Length - 6 ft (1.8 m) VALO Charger Power Supply is a Medical Grade Class II power supply and provides isolation from MAINS power	
VALO Charger	VALO 3.6VDC Lithium Iron Phosphate smart battery charger: <ul style="list-style-type: none"> <li>Automatic shut off when fully charged</li> <li>Auto-detection of defective batteries</li> <li>Protections: Thermal, Overcharge, Short-circuit, reverse polarity <ul style="list-style-type: none"> <li>Red LED – Charging</li> <li>Green LED – Empty or Fully Charged</li> <li>LED off – short circuit</li> </ul> </li> <li>Charging time: 1 – 3 hours</li> </ul> Rating: CE, WEEE		
VALO Batteries	Rechargeable: Safe chemistry Lithium Iron Phosphate (LiFePO4) RCR123A <ul style="list-style-type: none"> <li>Working Voltage: 3.2VDC, 400mAh (Watt hour needed)</li> </ul> Ratings: UL 62133-2, CE, RoHS, WEEE Cells are rated for >2000 charge/discharge cycles before hitting 70% of its capacity. It is advisable to replace batteries every 3–5 years, but life may vary depending on use. Replace as necessary.		
Operating, Charging Conditions	Temperature: +10°C to +32°C (+50°F to +90°F) Relative Humidity: 10% to 95% Ambient Pressure: 700 hPa to 1060 hPa		
Duty Cycle:	The curing light is designed for short-term operation. At maximum ambient temperature (32°C) 1 minute ON back-to-back cycling, 30 minutes OFF (cooling-off period).		

## Trouble Shooting

If the solutions suggested below do not rectify the problem, please call Ultradent at 800.552.5512. Outside the United States, call your Ultradent distributor or dental dealer.	
Problem	Possible Solutions
Light will not turn on	<ol style="list-style-type: none"> <li>Press the Time/Mode or Power Button to wake from Power Save Mode.</li> <li>Check the red Low Battery Indicator for battery charge status.</li> <li>Check that fresh batteries are correctly inserted into the unit.</li> <li>If red and yellow Warning LEDs are flashing this means the curing light has reached its internal temperature safety limit. Allow the curing light to cool down for 10 minutes or use a cool moist towel to cool the unit down quickly.</li> <li>If red Warning LED flashes and beeps continuously, call Ultradent Customer Service for repair.</li> </ol>
Light does not stay on for desired time	<ol style="list-style-type: none"> <li>Check Mode and Timing Lights for correct time input.</li> <li>Check the Low Battery Indicator for battery charge status.</li> <li>Check that fresh batteries are properly inserted into the unit.</li> </ol>
Light is not curing resins properly	<ol style="list-style-type: none"> <li>Check lens for residual cured resins/composites.</li> <li>Using proper amber UV eye protection, verify the LED lights are working.</li> <li>Check power level with light meter. If using a light meter, Ultradent recommends checking the curing light in Standard Power mode. NOTE: The true numeric output will be skewed due to the inaccuracy of common light meters and the custom LED pack the curing light uses. Light meters differ greatly and are designed for specific light guide tips and lenses.</li> <li>Check expiration date on curing resin.</li> <li>Ensure proper technique is being followed according to manufacturer recommendations.</li> </ol>

Batteries will not charge	<ol style="list-style-type: none"> <li>1) Make sure batteries are inserted into the charger in the correct orientation and allow batteries to charge for 1-3 hours.</li> <li>2) If red lights on the charger do not change to green, call Ultradent Customer Service to order replacement batteries and/or charger.</li> <li>3) If neither green nor red lights on the charger are visible, call Ultradent Customer Service to order or replace charger and/or AC adapter.</li> </ol>
Charger does not charge batteries	<ol style="list-style-type: none"> <li>1) Make sure charger is plugged in and AC adapter is plugged into a working power outlet.</li> <li>2) If green or red lights on the charger are not visible, call Ultradent Customer Service for new charger and/or AC adapter.</li> </ol>
Cannot change mode or time intervals	<ol style="list-style-type: none"> <li>1) Hold both Time/Mode and Power buttons down until a series of beeps indicates the curing light is unlocked.</li> </ol>

## 9. Miscellaneous Information

Guidance and Manufacturer's Declaration for Electromagnetic Emissions		
The curing light is intended for use in the electromagnetic environment specified below. The customer or user should ensure that it is used in such an environment.		
Emissions Test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The curing light uses batteries and is not affected by EMI, RF, or surge suppression. The curing light uses electrical and electromagnetic energy only for their internal functions. However, any RF emissions are very low and are not likely to cause interference in nearby electronic equipment.
RF emission CISPR 11	Class B	
Harmonic emissions IEC 61000-3-2	N/A	Harmonic emission and voltage fluctuation testing are not applicable to the curing light because it is battery powered. The curing light is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings for domestic use.
Voltage fluctuations/flicker emissions IEC 61000-3-3	N/A	

Guidance and Manufacturer's Declaration for Electromagnetic Immunity			
The curing light is intended for use in the electromagnetic environment specified below. The customer or user should ensure that it is used in such an environment.			
IMMUNITY test	IEC 60601 test level	Compliance level	Electromagnetic environment guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 8 kV contact ± 15 kV air	± 8 kV contact ± 15 kV air	Physical environment should be restricted to the following: 1. IP Code: IP20 2. Do not immerse in liquid. 3. Do not use around flammable gas. Unit is Non-APG and Non-AP. 4. Storage humidity range: 10% - 95% 5. Storage temperatures range: 10° C - 40° C
Electrical fast transient/ burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	± 2 kV for power supply lines Note 1: The curing light has no I/O ports	The curing light is battery powered and is not capable of connection to AC MAINS power.
Surge IEC 61000-4-5	± 1 kV line to line + 2 kV line to earth	± 1 kV line to line + 2 kV line to earth	Because the curing light is battery powered, it is not subject to electric transients, surges, voltage dips, shorts, interruptions, or variations on AC MAINS power.
Voltage, dips, shorts, interruptions, and variations on the power supply input lines  IEC 61000-4-11	<5% U (>95% dip in U for 0.5 cycle)  40% U (60% dip in U for 5 cycles)  70% U (30% dip in U for 25 cycles)  <5% U (>95% dip in U for 5 s)	<5% U (>95% dip in U for 0.5 cycle)  40% U (60% dip in U for 5 cycles)  70% U (30% dip in U for 25 cycles)  <5% U (>95% dip in U for 5 s) Note 2: Self recovers	The accessory battery charger may be subject to the above but it is separate from and not critical to the operation of the VALO Cordless curing light.  If battery voltage of the curing light drops to 4VDC the unit will not allow operation. The VALO Cordless curing light will turn off. When new batteries are put in and the proper power levels are restored, the curing light will restart and return to the same state before power loss. The curing light will self-recover in the event of power loss.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical, residential, home health care, commercial, hospital, or military environment.
Note 1: The curing light is not equipped with any ports or any accessible I/O lines.			

Guidance and Manufacturer's Declaration for Electromagnetic Immunity for non-life support systems			
The curing light is intended for use in the electromagnetic environment specified below. The customer or user should ensure that it is used in such an environment.			
IMMUNITY test	IEC 60601 test level	Compliance level	Electromagnetic environment guidance
Conduction RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3 Vrms 150 kHz to 80 MHz	Portable and mobile RF communications equipment should be used no closer to any part of the curing light, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance: $d = \left[ \frac{3,5}{V_1} \right] \sqrt{P}$ $d = \left[ \frac{3,5}{E_1} \right] \sqrt{P} \quad 80 \text{ MHz to } 800 \text{ MHz}$ $d = \left[ \frac{7}{E_1} \right] \sqrt{P} \quad 800 \text{ MHz to } 2,5 \text{ GHz}$ P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol: 
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m 80 MHz to 2.5 GHz	

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

a. Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the curing light is used exceeds the applicable RF compliance level above, the VALO Cordless curing light should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the curing light.

b. Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

**Guidance and Manufacturer's Declaration for recommended separation distances between portable and mobile RF communications equipment and the VALO Cordless curing light**

The curing light is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The user of the curing light can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the curing light as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter (P in Watts)	Separation distance according to frequency of transmitter (meters)		
	150 kHz – 80 MHz $d = \left[ \frac{3.5}{V_1} \right] \sqrt{P}$	80 MHz – 800 MHz $d = \left[ \frac{3.5}{E_1} \right] \sqrt{P}$	800 MHz – 2.5 GHz $d = \left[ \frac{7}{E_1} \right] \sqrt{P}$
0.01	0.12 meters	0.035 meters	0.07 meters
0.1	0.37 meters	0.11 meters	0.22 meters
1	1.7 meters	0.35 meters	0.7 meters
10	3.7 meters	1.11 meters	2.22 meters
100	11.7 meters	3.5 meters	7.0 meters

The curing light has been tested according to IEC 60601-1-2:2014 and passed under radiated field strengths of 10 V/m between 80-MHz to 2.5 GHz. The value of 3Vrms corresponds to V1 and the value 10 V/m corresponds to E1 in the formulas above.

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Report any serious incident to the manufacturer and the competent authority.

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
MEDICAL - GENERAL MEDICAL EQUIPMENT  
AS TO ELECTRICAL SHOCK, FIRE, AND  
MECHANICAL HAZARDS ONLY  
IN ACCORDANCE WITH  
ANSI/AAMI ES 60601-1 (2005) + AMD (2012)  
CAN/CSA-C22.2 NO. 60601-1 (2008) +  
(2014)  
E357088

EN - For professional use only  
DE - Nur zur Anwendung durch den Zahnarzt  
FR - Pour usage professionnel seulement  
NL - Alleen voor professioneel gebruik  
IT - Solo per uso professionale  
ES - Solo para uso profesional  
PT - Apenas para utilização profissional  
SV - Endast för professionell användning  
DA - Kun til professionel brug  
FI - Vain ammattikäyttöön  
EL - Για επαγγελματική χρήση μόνο  
CS - Pouze pro odborné použití.  
BG - Само за професионално ползване  
HR - Samo za profesionalnu uporabu  
ET - Ainult professionaalseks kasutamiseks  
HU - Kizárólag professzionális felhasználásra  
GA - Le haghaidh úsáide gairmíola amháin  
LV - Tikai profesionālai lietošanai  
LT - Tik profesionaliam naudojimui  
MT - Għall-użu professjonali biss  
RO - Numai pentru uz profesional  
SK - Len na profesionálne použitie  
SL - Samo za profesionalno uporabo  
NO - Kun for profesjonell bruk  
PL - Wyłącznie do użytku profesjonalnego.  
RU - Только для профессионального использования  
TR - Sadece profesyonel kullanım için  
AR - لايسخدم من قبل متخصصين فقط  
TH - สำหรับการใช้งานโดยมืออาชีพเท่านั้น  
KO - 전문가 전용  
ZH - 仅供专业人士使用

EN - Importer  
DE - Importeur  
FR - Importateur  
NL - Importeur  
IT - Importatore  
ES - Importador  
PT - Importador  
SV - Importör  
DA - Importør  
FI - Maahantuoja  
EL - Εισαγωγέας  
CS - Dovezce  
BG - Вносител  
HR - Uvoznik  
ET - Importija  
HU - Importőr  
GA - Allmhairéoir  
LV - Importētājs  
LT - Importuotojas  
MT - Importatur  
RO - Importator  
SK - Dovezca  
SL - Uvoznik  
NO - Importør  
PL - Importer  
RU - Импортёр  
TR - İthalatçı  
AR - الممور  
TH - ผู้นำเข้า  
KO - 수입자  
ZH - 进口商

EN - Medical Device  
DE - Medizinprodukt  
FR - Dispositif médical  
NL - Medisch instrument  
IT - Dispositivo medico  
ES - Dispositivo médico  
PT - Dispositivo médico  
SV - Medicinska naprava  
DA - Mediskinsk anordning  
FI - Lääketieteellinen laite  
EL - Ιατροτεχνολογικό προϊόν  
CS - Lékařské zařízení  
BG - Медицинско изделие  
HR - Medicinski uređaj  
ET - Meditsiinsiseade  
HU - Orvostechnikai eszköz  
GA - Uirlis leighis  
LV - Medicīniskā ierīce  
LT - Medicinos prietaisai  
MT - Apparat mediku  
RO - Dispozitiv medical  
SK - Lekárský prístroj  
SL - Medicinska naprava  
NO - Medisinsk enhet  
PL - Wyrób medyczny  
RU - Медицинское изделие  
TR - Tıbbi Cihaz  
AR - جهاز طبي  
TH - อุปกรณ์ทางการแพทย์  
KO - 의료 기기  
ZH - 医疗器械

[www.ultradent.com](http://www.ultradent.com) / 1-800-552-5512 / 801-572-4200  
Report any serious incident to the manufacturer and the competent authority.

© 2024 Ultradent Products Inc.

---

Manufactured by

 Ultradent Products, Inc.  
505 West Ultradent Drive (10200 South)  
South Jordan, UT 84095  
88950AR13 112024

 **ULTRADENT**  
PRODUCTS, INC.

 Ultradent Products GmbH  
Am Westhoyer Berg 30  
51149 Cologne Germany

Australian Sponsor:  
Ultradent Australia Pty Ltd  
22/2 Market St. Sydney, NSW 2000  
1800 29 09 29 | [www.ultradent.com.au](http://www.ultradent.com.au)