CE

C-Soft 2.6/C-Soft 2.8

4X180=V2.6/ 4X180=V2.8

[DE]	Gebrauchsanweisung (Fachpersonal)	3
[EN]	Instructions for use (Qualified Personnel)	25
[FR]	Instructions d'utilisation (Personnel spécialisé)	47
[IT]	Istruzioni per l'uso (Personale tecnico specializzato)	71
[ES]	Instrucciones de uso (Personal técnico especializado)	93
[PT]	Manual de utilização (Pessoal técnico)	115
[NL]	Gebruiksaanwijzing (Vakmensen)	137
[SV]	Bruksanvisning (Fackpersonal)	159
[JA]	取扱説明書 (有資格担当者)	181

Inhalt

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Table of contents

1 Introduction Ottobock C-Soft.....	25
1.1 Intended Use	25
1.2 Qualification of the Prosthetist	25
1.3 Safety Instructions	25
2 Scope of Delivery and Accessories	26
3 Installation and Pairing	27
3.1 Hardware requirements	27
3.2 Software requirements	27
3.3 Installation	28
3.4 Uninstalling	30
3.5 Pairing 60X4 BionicLink PC	31
3.6 Starting the Program	31
3.7 Closing the program.....	32
3.8 Further information.....	32
4 Using the Ottobock Data Station.....	32
4.1 General Description of the Program's User Interface	32
4.2 First Steps with the Ottobock Data Station	33
4.3 Job Selection	34
5 Using Ottobock C-Soft.....	35
5.1 General Description of Ottobock C-Soft	35
5.2 Selecting the Type of Connection and Establishing the Connection.....	36
5.3 Adjustments on the Knee Joint	39
5.4 Data Overview.....	41
6 Liability Exclusion.....	42
7 Ottobock C-Soft Licence Terms.....	42
8 Symbols on the product.....	42
9 Ottobock C-Soft Licence Terms.....	43

INFORMATION**Last update:** 2016-03-10

- Please read this document carefully before using the product.
- Follow the safety instructions to avoid injuries and damage to the product.
- Please contact the manufacturer if you have questions about the product (e.g. regarding the start-up, use, maintenance, unexpected operating behaviour or circumstances). Please see the back for contact information.
- Please keep this document in a safe place.

1 Introduction Ottobock C-Soft

Ottobock has advanced its software for adjusting the electronic leg prosthesis systems (C-Leg, C-Legcompact, Ottobock Compact). With C-Soft, prosthetists can make use of state-of-the-art application technology for individual adjustments of the electronic leg prosthesis systems.

The Ottobock C-Soft program offers prosthetists the possibility to optimally adjust the electronic leg prosthesis system to a patient. It guides users step by step through the adjustment procedure and assists them with images and text. After the adjustment, users can save the setting and make a print-out for documentation purposes. If need be, these values can be loaded and fed into the electronic leg prosthesis system.

1.1 Indication for use

The Ottobock C-Soft adjustment software 4X180=V2.6 /4X180=V2.8 is used for adjusting the parameters of electronic Ottobock knee joint systems with the following article numbers: 3C95/3C85, 3C96/3C86, 3C96-1/3C86-1, 3C98/3C88, 3C98-1/ 3C88-1, 3C98-2/3C88-2.

1.2 Qualification of the Prosthetist

The adjustment software may only be used by users who have participated in the relevant product training and have been certified for the application. Additional product training courses may become necessary to qualify for software updates.

1.3 Safety Instructions

Explanation of Symbols

△ CAUTION Warnings regarding possible risks of accident or injury.**NOTICE** Warnings regarding possible technical damage.**△ CAUTION**

Alignment and adjustment error. During alignment and adjustment of the prosthesis, errors can occur resulting in malfunction of the joint up to loss of function due to structural failure. This can cause the patient to fall.

- Participation in an Ottobock product training course for the C-Leg / C-Leg compact / Ottobock Compact is obligatory before the first fitting. Additional product training courses may become necessary to qualify for fitting product updates.
- Correct maximum toe load setting is an important criterion for the quality of the fitting. If the maximum toe load has been adjusted too high, it may happen that the prosthesis does not switch to the swing phase. If the maximum toe load has been adjusted too low, it may happen that the prosthesis switches to the swing phase too early.
- During the data transfer (PC to knee joint), the patient must remain standing or sitting still and the BionicLink or the communication cable must not be removed.
- If the prosthesis wearer uses crutches or walking canes during adjustment, readjustment is required as soon as he stops using these walking aids.

△ CAUTION

Risk of falling while making adjustments. If the connection to C-Soft is interrupted or established unintentionally while making adjustments, this can lead to unexpected changes of movements. This can cause the patient to fall.

- The knee joint must always be disconnected from C-Soft after adjustments have been completed.
- As long as the patient is connected to Ottobock C-Soft when making adjustments, he or she must not be left unattended.
- If the connection fails while making adjustments, the prosthetist must immediately warn and secure the patient.

NOTICE

Damaged contacts. The elastomer cap serves to protect the electrical contacts of the C-Legcompact/Ottobock Compact. Attaching the elastomer cap improperly can damage the electrical contacts.

Carefully put on the elastomer cap after making adjustments (only with C-Leg compact 3C96/3C86 and Ottobock Compact 3C95/3C85).

INFORMATION

Familiarization with the prosthesis. After the amputee has got accustomed to the prosthesis, his gait pattern will change. For this reason it is necessary to repeat the complete adjustment procedure approximately one week after the first fitting. The process of habituation, however, will take longer than one week. Therefore, it is to be expected that the patient's activity will change.

2 Scope of Delivery and Accessories

- 1 pc. 4X180=V2.6 C-Soft 2.6 software CD
- Update to 4X180=V2.8 C-Soft 2.8 via the Ottobock support website (see section 3.8.2 Software updates)
- Instructions for Use 647G268

C-Soft supports the following PC interface adapters:

- 60X1 USB-LWL (fibre optic cable) Adapter
- 60X2 RS232-LWL (fibre optic cable) Adapter
- 60X4 BionicLink PC Adapter with 60X3 Bionic Link (for Windows XP and Windows Vista only)
- 60X5 BionicLink PC Adapter with 60X3 Bionic Link

INFORMATION For explanations on the installation of the 60X1 USB-LWL Adapter, 60X2 RS232-LWL Adapter, and 60X5 BionicLink PC Adapter, please refer to the product-specific instructions for use.

Product	Instructions for Use
60X1	647H481
60X2	647H496
60X3	647G192
60X5	647G506

3 Installation and Pairing

3.1 Hardware requirements

3.1.1 Recommended PC hardware requirements

For 32-bit platforms:

- PC with a Pentium IV 650 / 3.4 GHz processor (x86) (or faster)
- Minimum of 1 GB RAM (working memory)

For 64-bit platforms:

- PC with an Intel Core2 (x64) or faster processor
- Minimum of 2 GB RAM (working memory)

For 32 bit and 64 bit:

- Minimum 1 GB available hard disk space
- Graphics card with Open GL support
- Minimum resolution 1024 x 768 at 96 DPI (higher DPI may cause display errors, depending on the resolution), 32 bit colour depth (16.7 million colours)
- CD-ROM or DVD-ROM drive
- 1 available USB connection
- Mouse and keyboard
- Internet connection

3.2.1 Minimum PC hardware requirements

The minimum hardware requirements depend on the minimum hardware requirements of the operating system. If the following requirements are greater, then these apply:

- PC with a Pentium III / 1 GHz 32-bit (x86) or faster processor
- 512 MB RAM (working memory)
- 1 GB available hard disk space
- Graphics card with Open GL support
- Resolution 1024 x 768 at 96 DPI, 32 Bit colour depth (16.7 million colours)
- CD-ROM drive
- 1 available USB connection
- Mouse and keyboard

3.2 Software requirements

The following 32 bit operating systems are supported:

- Microsoft Windows XP Professional, SP2 or higher (recommended: SP3)
- Microsoft Windows XP Home Edition, SP2 or higher (recommended: SP3)
- Microsoft Vista, all editions except Microsoft Vista Starter (recommended: SP2)
- Microsoft Windows 7, all editions except Microsoft Windows 7 Starter
- Microsoft Windows 10

The following 64-bit operating systems are supported:

- Microsoft Windows 7, all editions except Microsoft Windows 7 Starter
- Microsoft Windows 8
- Microsoft Windows 10

INFORMATION

Please note that the latest USB driver version must always be installed for the Bluetooth connection via the 60X5 BionicLink PC. You can download the latest driver version from the following site: <http://ottobock.com/en/datasheet>. The latest version can also be downloaded via the "software update" function (3.8.2 Software Updates).

3.2.1 Additional recommended software components:

The following additional software components are recommended:

- Internet browser
- A MAPI-compatible e-mail system (e.g. Microsoft Outlook, Mozilla Thunderbird ...)

3.2.2 Additional installed software components:

The following additional software components are required and are installed during installation of the software (if they are not installed on the computer already):

- Microsoft XML 6.1
- Microsoft .NET Framework 2.0

3.3 Installation

To switch the user interface language to Turkish, the C-Soft software update to Version V2.8 must first be downloaded and installed. For information on downloading the update, see chapter "Software updates".

INFORMATION

Microsoft Windows software user rights

- Administrator rights are required to install the software.
- At least standard user rights are required to use the software.

To run the installation program, start Windows and insert the Ottobock C-Soft 2.6 CD into the CD-ROM drive.

Installation by CD is started automatically after the CD has been inserted.

INFORMATION

If the installation does not start after the CD is inserted

Windows XP: open "My computer" (key combination [Windows logo key + E]). Click on the CD drive in the left window and start the program "setup.exe" in the right window by double-clicking on it.

Windows 7/Windows 8: open "Computer" (key combination [Windows logo key + E]). Right-click on the CD drive in the right window. Click on the menu item "Install or run program from your media" with the left mouse button.

Windows 10: open "Computer" (key combination [Windows logo key + E]). Right-click on the CD drive in the left window. Click on the menu item "Install or run program from your media" with the left mouse button.

After starting the setup program, the Setup Wizard will take you through a number of screens asking you to:

- 1) Read and accept the licence agreement.

2) Enter your user information and click "Continue".

3) You can register the software with Ottobock. In order to do so, enter your company information (fields marked with an asterisk* must be completed) and click "Continue".

If you do not want to register, click "Skip registration".

4a) **To register:** After entering your company information, click "Continue". You will then receive a summary of your data that will be transmitted to Ottobock during the registration process.

- Click "E-mail registration". Your data will be e-mailed to Ottobock.
- If you do not have access to e-mail but do have internet access on your PC, click the "Copy to clipboard" button. Afterwards, email the contents of the clipboard by means of a webmail service (e.g. T-Online, Yahoo, ...) to **registration.austria@ottobock.com**.
- Alternatively, you receive a printout you can fax to Ottobock when you click the "Print" button.
Please send this fax to: +43-1 526 79 85

Otto Bock Healthcare Products GmbH - Kaiserstraße 39 - 1070 Wien - Austria

- 4b) **To skip registration:** After clicking "Skip registration", the "Ready to install the programme" window is displayed.

- Ottobock C-Soft 2.6/C-Soft 2.8 can be accessed after having started the Ottobock Data Station.
- **[INFORMATION]** The Ottobock Data Station is a superior software platform that is used to start the respective programs. The installation program therefore does not create an icon for Ottobock C-Soft 2.6/C-Soft 2.8.

3.4 Uninstalling

Windows XP

1. Open the Control Panel:

Click "Start > Control Panel > Software".

→ All installed software packages are displayed.

2. Select the entry "Ottobock C-Soft".

3. Click the "Remove" button.

→ All "Ottobock C-Soft" software components will be removed from the hard drive.

However, the entered patient data and the Data Station are not deleted.

Windows 7

1. Open the Control Panel:

Click "Start > Control Panel > Programs and Features".

→ All installed software packages are displayed.

2. Select the entry "Ottobock C-Soft".

3. Click the "Uninstall" button.

→ All "Ottobock C-Soft" software components will be removed from the hard drive.

However, the entered patient data and the Data Station are not deleted.

- 5) By clicking the "Install" button, the installation process is started. The installation can take a few minutes.

- 6) Afterwards, a window appears informing you about the end of the installation process.

- 7) Click the "Finish" button to complete the installation. An icon for starting the Ottobock Data Station has been created on the Windows desktop.

Windows 8/Windows 10


1. Press the key combination [Windows logo key  +X]
2. Click the "Programs and Features" entry.
→ All installed software packages are displayed.
3. Select the entry "Ottobock C-Soft".
4. Click the "Uninstall" button.
→ All "Ottobock C-Soft" software components will be removed from the hard drive.
However, the entered patient data and the Data Station are not deleted.

3.5 Pairing 60X4 BionicLink PC

INFORMATION

Pairing, i. e. first-time connection of the **PC** to the Bluetooth adapter, is only **required** with the **60X4 BionicLink PC**. With the **60X5 BionicLink PC**, the connection selection is made when establishing the connection in C-Soft.
BionicLink PC **60X4** is not compatible with Windows 7/Windows 8/Windows 10 .

1. In the operating system, open "*Bluetooth Places*".
2. In "*Bluetooth Places*" search for available Bluetooth Services.
3. Before selecting a Bluetooth™ Service, check the device number on the back side of the housing of the **60X3 BionicLink**. Select the Bluetooth™ Service "*BionicLink XXX*" with the same device number.

4. Double-click the "*SPP slave on BionicLink XXX*" symbol to establish the connection.
5. The operating system informs you that the "*Bluetooth PIN code is required*".
6. Enter the Bluetooth **PIN-Code "2534"**.
7. In case of successful connection establishment, a dialogue window appears displaying the device number and the **COM-PORT** number. The **COM-PORT** number should be kept in mind, as you will need it during the establishment of the connection in C-Soft for the **COM-PORT** selection.

3.6 Starting the Program

The installation program adds an entry in the Windows Start menu and on the Windows desktop. The suggested standard path is *Program files* → *Ottobock* → *Ottobock Data Station*.

3.7 Closing the program

You can close Ottobock C-Soft using one of four methods:

- "Exit" command on the "File" menu
- "Close Window" button
- Selecting the Windows system menu "Close via Start Menu"
- Using the key combination Alt + F4.

3.8 Further information

3.8.1 Obtaining Program Version Information

The current program version of the installed software can be obtained from the menu bar. You should have this information available whenever you contact the software support.

1. Click "Help > About" in the menu bar.
→ A window with the currently installed software will open.
2. Click on the desired entry to obtain further information.


3.8.2 Software updates

Software and driver updates will be made available for free at regular intervals.

1. Click "Help > Otto Bock Support Website" in the menu bar when you are connected to the Internet.
→ The webpage for the "Data Station" software will open.
2. Search for the description of the software that is to be updated under "All".
3. Search for the applicable version.
4. Click "Download" in the column on the right in order to download the update.
5. Extract the "ZIP file" and execute it.

4 Using the Ottobock Data Station

4.1 General Description of the Program's User Interface

The start tab of the Ottobock Data Station consists of the following two sections:

- Overview of companies and patients (group list and patient list) **A** ;
- Overview of patient data (detail list) **B** .

4.2 First Steps with the Ottobock Data Station

1. Changing the language in the software (optional):

In the menu "Tools", select "Settings". In the next window, select the tab "Customization". In the pull-down menu, select the desired language and click the button "OK" to confirm.

2. Selecting the company

Select the "company" by clicking with the left mouse button on the company icon or on the "company" name (highlighted item **1**).

3. Selecting or creating a patient file

- To select an **existing patient file**, click on the patient name (highlighted item **2**).

The patient data overview displays all personal data that have already been entered (and which can be changed, if need be) as well as the jobs performed.

- To **create a new patient file**, click the button "New patient" (highlighted item **3**).

In the patient data overview, enter the personal data of the patient (last name, first name, etc.).

*Fields marked with ** are mandatory fields.*

4. Editing or creating a new job

- To **create a job**, click the button "New job" (highlighted item **5.1**).

Select "C-Soft" as type of job.

- To **edit an existing job**, double-click the job name in the list of "Jobs" (highlighted item **4**) (The selected job must have been initially edited in C-Soft). C-Soft will start automatically.

5. Changing to job selection (only when creating a new job)

Select the tab "New job" to get to the program selection (highlighted item **5.2**).

6. Saving and printing the settings

Click the save symbol (highlighted item 1) and save the job under the desired name.

Click the print symbol (highlighted item 2) and select the printer.

Saving and printing is possible at any time.

4.3 Job Selection

Click the "C-Soft" button (highlighted A) in the lower extremity selection section to start C-Soft.

Entering the password for the first time

Password (Unlock-PIN) entry is required after "Ottobock Data Station – C-Soft" has been started. If the software is used for the first time, the password (Unlock-PIN) you received within the scope of your successful certification must be entered. If you are certified in C-Soft V1.0/V1.2/V2.0/V2.2/V2.4, the password (unlock pin) you received at that time is also valid for version 2.6 and version 2.8.

After entry of the password (Unlock-PIN), you will be requested to enter and confirm another password (Unlock-PIN) which this time can be freely chosen by you.

This application requires a password. Please enter a valid username and password.

INFORMATION The password (Unlock-PIN) you received within the scope of certification can be used at any time instead of your personal password. If this has been done, you will again be requested to enter and confirm another password which can be freely chosen by you.

5 Using Ottobock C-Soft

5.1 General Description of Ottobock C-Soft

Navigation section (highlighted item A):


All adjustments can be made step by step, one after the other. Parameters that cannot be adjusted are hidden. Click the different tabs to change between the adjustment steps. Very experienced users can also use the expert mode ("Advanced Settings").

C-Soft can be closed at any time by clicking the button (top right). The start screen of the Data Station opens automatically.

Adjustment section (highlighted item B):

All adjustment steps are explained with text and images. Follow the software instructions for adjusting the joint.


Status section (highlighted item C):

The status of the connection, battery, as well as knee angle, maximum toe load, dynamic toe load and heel load are visualized in real time. To each adjustment step, additional useful information can be displayed by clicking the "info" button (highlighted item 1). The battery symbol (highlighted item 2) indicates the

remaining battery capacity of the knee joint. The connection symbols indicate the current connection status. If the symbol is red (highlighted item 3.1), there is no connection between the adjustment software and knee joint. If the symbol is green (highlighted item 3.2), a connection has been established and adjustments can be made.

The presentation of the knee joint visualizes the current knee angle in real-time (highlighted item 4).

The presentation of the prosthetic foot visualizes the currently set maximum toe load value (highlighted item 5), the current toe load value (highlighted item 6) and the current heel load value (highlighted item 7) in real-time.

When being in the status section, it is also possible to change to the next adjustment step (by clicking on the respective button).

5.2 Selecting the Type of Connection and Establishing the Connection

Selection of the type of connection depends on the adapter used. The different adapters must be installed before starting the program.


Section 2 "Scope of Delivery and Accessories" contains an overview of the interface adapters.

- **Ottobock Bluetooth Adapter (*) (highlighted item 1)**:

When using a wireless Bluetooth adapter, select the tab "Ottobock Bluetooth Adapter" (highlighted item 1).

With 60X5 BionicLink PC:

1. Click the button "Connect via Bluetooth" (highlighted item 2).
2. In the Bluetooth device selection window, select the desired BionicLink device according to the number on the back side of the housing (highlighted item 3), and click the "OK" button (highlighted item 4) to confirm. If only one BionicLink device is available, the selection will automatically be skipped.
3. In case of successful connection establishment, the connection symbol will change to "connected".

Reconnect:

Click the button "Reconnect" (highlighted item 5) to reopen the Bluetooth device selection window and select another BionicLink. If only one BionicLink device is available, the selection will automatically be skipped (see item 1).


* Wireless Bluetooth adapters are:


- 60X3 BionicLink with 60X5 BionicLink PC
- 60X3 BionicLink with 60X4 BionicLink PC

With 60X4 BionicLink PC:

1. Click on "Optional: Manual COM-Port selection" and select the COM-Port that had been displayed during "pairing" of the 60X4 BionicLink PC (see section 3.6) (highlighted item 2.2).
2. Click the button "Connect via Bluetooth" (highlighted item 2.1). In case of successful connection establishment, the connection symbol will change to "connected".

Connection failed::

Measures:

1. Verify whether all plug connections have been established:
 - USB plug on the PC
 - BionicLink plug on the knee joint
2. Check the status LEDs of the adapters (refer to the instructions for use of the respective component).
3. With 60X4 BionicLink PC, "pairing" must be carried out before the initial operation (see section 3.6).
4. Verify, whether 60X3 BionicLink is within reach (10 metres).



With Ottobock USB/Serial Adapter (**)(highlighted item 1.2):

When using a serial adapter, select the tab "Ottobock USB/Serial Adapter" (highlighted item 1.2). The installed and operational serial interface adapter is displayed.


Click the button "Connect" (highlighted item 2.1). In case of successful connection establishment, the connection symbol will change to "connected".

** Serial adapters are:

- 60X1 USB-LWL (fibre optic cable) Adapter
- 60X2 RS232 Adapter

Connection failed:

Measures:

1. Verify whether all plug connections have been established:
 - USB/RS232 plug on the PC
 - LWL (fibre optic cable) plug on the knee joint
2. Check the status LEDs of the adapters (refer to the instructions for use of the respective component).

To get to the next step:

After successful connection establishment, click the "Calibration" tab in the navigation section to change to calibration.

5.3 Adjustments on the Knee Joint

INFORMATION

The program does not effect any automatic changes to the settings!
Each change must be made or confirmed by the prosthetist!

C-Soft describes the proceeding for each adjustment step in the adjustment section. Additional information is available in the status section ("info" button).

Calibration (Zero Setting)

INFORMATION

Instruct the patient in accordance with the explanations in the software's adjustment section.

Calibration

The amputee is standing between parallel bars or on crutches with the prosthesis fully extended and the prosthetic foot unloaded (made possible by moving the residual limb backwards).
The prosthetic foot should not touch the ground.

The sensors of the prosthesis have to be calibrated every time a connection is established.
To do so, press the button "Calibrate".
After successful calibration you will hear several beeps from the knee joint.

Knee Angle

Toe- and heel load

Calibrate

Adjustment 2

Calibration is required for all electronic knee joints compatible with C-Soft. **Adjustments are not possible without prior calibration.**


Click the "Calibrate" button (highlighted item 1) to calibrate the knee joint.

INFORMATION Calibration is required each time a connection has been established.

Display in the status section:

Zero setting OK:

All functions of the program can be used.

Zero setting not OK:

The knee angle has not reached its zero value.

Possible causes:

Wrong patient position or mechanical blockage due to: caught fabric, caught cosmetic foam cover, etc.

Zero setting not OK:

It was not possible to calibrate the tube with strain gauge.

Possible causes:

The tube with strain gauge was not connected or was not calibrated correctly, or the foot was loaded.

To get to the next step:

After successful calibration, click the "Adjustment" tab in the navigation section or the "Adjustment" button (highlighted item ②) to change to the knee joint settings.

All adjustments are explained in the adjustment section as well as in the navigation section ("Info" button) in C-Soft (see highlighted items in Section 5.1).

5.4 Data Overview

The data overview displays the currently adjusted values and the serial number. If values were already saved for a patient, they will be displayed next to the currently adjusted values. The currently adjusted values and the values already saved can be distinguished by the following symbols.

Values currently adjusted in the electronic knee joint:

Values saved:

Save (Load) data

Data is backed up on your computer and can be loaded at any time.

By clicking on 'Save current values', you will be prompted to enter a job name. Enter the job name (e.g. date or trial fitting) and confirm.

1. Saving (currently adjusted values)

To save the values currently adjusted in the electronic knee joint, click on the field 'Save current values'.

Save values currently adjusted in the electronic knee joint:

2. Loading (data already saved)

To load saved values, click on the field 'Load saved values'.

Load saved values into the electronic knee joint:

6 Trademarks

All product names mentioned in this document are subject without restriction to the respective applicable trademark laws and are the property of the respective owners.

All brands, trade names or company names may be registered trademarks and are the property of the respective owners.

Should trademarks used in this document fail to be explicitly identified as such, this does not justify the conclusion that the denotation in question is free of third-party rights.

- Microsoft, Outlook and Windows are registered trademarks of Microsoft Corporation.
- Pentium, Intel and Intel Core are registered trademarks of Intel Corporation.

7 Liability Exclusion

The manufacturer explicitly states that this product may only be used in combination with components that were authorized by the manufacturer (see instructions for use and catalogs). The manufacturer does not assume liability for damage caused by component combinations which it did not authorize.

8 Symbols on the product

Legal manufacturer

Declaration of conformity according to the applicable European directives

9 Ottobock C-Soft Licence Terms

Manufacturer: Otto Bock Healthcare Products

PLEASE READ THESE END USER LICENCE AGREEMENT ("EULA" OR "LICENCE AGREEMENT") CAREFULLY BEFORE INSTALLING OR PUTTING THE OTTO BOCK SOFTWARE INTO OPERATION. BY INSTALLING OR PUTTING THE OTTO BOCK SOFTWARE INTO OPERATION, YOU DECLARE YOUR AGREEMENT WITH THE TERMS OF THIS EULA.

IF YOU HAVE NOT UNDERSTOOD THE TERMS OF THIS LICENCE AGREEMENT, DO NOT USE THE OTTO BOCK SOFTWARE. IF YOU DO NOT AGREE TO THE TERMS OF THIS LICENCE AGREEMENT, YOU CAN GIVE THE OTTO BOCK SOFTWARE BACK TO WHERE YOU HAVE THE SOFTWARE ACQUIRED FOR A REFUND OF THE PURCHASE PRICE. IF YOU HAVE ACQUIRED THE OTTO BOCK SOFTWARE COMBINED WITH A HARDWARE DEVICE, YOU CAN GIVE BACK BOTH THE HARDWARE AND THE SOFTWARE FOR A REFUND OF THE PURCHASE PRICE.

IF YOU DO NOT AGREE TO THE TERMS OF THIS EULA, YOU ARE NOT AUTHORIZED TO INSTALL AND/OR USE THE OTTO BOCK SOFTWARE.

1. Parties and Subject Matter

- 1.1 By accepting the EULA, an Agreement is made between the above Ottobock Company and you; if you install or use the Ottobock software for the purpose of a company (e.g. as an employee of the company), then the Agreement is made with the company (so that "you" refers to the company). In case of doubt, it shall apply that the Agreement has been made with you.
- 1.2 The licence shall comprise the above Ottobock software inclusive of related documentation (hereinafter referred to as "Ottobock Software"). The terms of this EULA shall also apply to any Ottobock Software updates replacing and/or complementing the original software unless such an update is subject to a separate licence agreement. Any maintenance of the Ottobock Software or other Ottobock products is not subject matter of this Agreement, but shall be covered by a separate agreement if need be.
- 1.3 Insofar as the Ottobock Software is to be used for providing guideline proposals for patient's fitting (e.g. for special adjustments of Ottobock products according to measurements taken by the customer), it should be considered that these proposals are Ottobock's recommendations without obligations. Your own professional responsibility for the proper patient's fitting remains unaffected. Therefore, Ottobock will do its best to provide proposals and correct recommendations, but cannot be held liable for them. Accordingly, Ottobock cannot be held liable even for possible errors in the recommendations made, except for intent.
- 1.4 Insofar, when the Ottobock Software trainings are required for the use, the software may only be used by the appropriately trained personnel.

2. Licence

- 2.1 The Ottobock Software is protected in many countries by national copyright laws and international agreements. The cases of breach of these proprietary rights will be pursued by Ottobock. You are entitled to use the Ottobock Software provided that it was stipulated by the law or these rights were expressly granted to you according to this Licence Agreement. Ottobock and/or the licensor(s) of Ottobock remain owners of all proprietary rights or other rights to the Ottobock Software.
- 2.2 By accepting this EULA, you are granted non-exclusive right to use the Ottobock Software as intended and for own needs. Own needs mean needs of the particular contract partner only. In case of medical supply stores or other companies providing services to patients related to Ottobock products (e.g. adjustment of prostheses), the own needs in terms of this licence mean also needs of the particular patient who wants to make the proper use of the Ottobock products.
- 2.3 You are authorised to install the Ottobock software on several computers at the same time, or to use it on a network where the software can be accessed on more than one computer simultaneously. Making the software accessible to the public is not permitted. In terms of this EULA, Ottobock products that use the software in accordance with this agreement are also computers.
- 2.4 You are not authorized to modify, process, translate, arrange or adapt in any manner whatsoever and/or decompile, disassemble or subject the Ottobock Software to reverse engineering unless expressly allowed to do so by the law.

3. Sublicensing, Transfer

Granting of sublicenses for using the Ottobock Software is not permitted. Renting, leasing, loaning, or otherwise temporarily ceding the Ottobock software to a third party is also prohibited. A permanent transfer of rights to use you were granted is only possible provided that:

- the third person to which you transfer the rights has accepted the terms of this Licence Agreement;
- the Ottobock Software will be transferred as a whole and not only its parts, and
- you do not keep any copies of the Ottobock Software, even not in random access or cache memories of your computer.

4. Warranty

A warranty for the Ottobock Software complies with the terms based on which the Ottobock Software was ceded to you. As a rule, these comprise the applicable General Terms of Business of Ottobock as well as legal provisions. This Licence Agreement does not construe the right to raise additional claims against Ottobock in case of defects.

5. Limitation of Liability

Ottobock is liable to you for a culpable breach of so called cardinal obligations to reimbursement of expenses and damage compensation (referred to in the following item 5 "damage compensation") in accordance with the legal regulations. Cardinal obligations mean all obligations, the breach of which endangers the accomplishment of the purpose of the agreement, as well as all obligations, which may only be met provided that the Agreement is performed properly and that you can regularly rely on their observation. However, if the breach of cardinal obligations results from slight negligence only and does not result in injury of life, body or health, your claims for damage compensation will be limited to the typically foreseeable damage.

Additionally, Ottobock will be liable in accordance with the German Product Liability Act; in cases of intention and gross negligence; for injury of life, body or health; in case of a warranty by Ottobock as well as all other cases of imperative liability in accordance with the legal regulations.

For the rest, your claims for damage compensation against Ottobock, regardless of the legal basis, especially due to material defect, legal defect and/or breach of other contractual obligations by Ottobock, its legal representatives, employees or persons acting on its behalf or due to unauthorized handling, are excluded.

Insofar as the Ottobock's liability is limited or excluded according to the foregoing provisions, the same applies to the personal liability of the legal representatives, employees or persons acting on Ottobock's behalf.

The foregoing provisions do not imply a change in the burden of proof to your disadvantage.

Your claims for damage compensation are subject to statute of limitations periods and their commencement as specified by law.

6. Governing Law and Jurisdiction

This Licence Agreement shall be subject to the manufacturer's law. The United Nations Convention on Contracts for International Sale of Goods does not apply to this Agreement. Exclusive place of jurisdiction for all disputes resulting from or in relationship to this Licence Agreement is the location of the manufacturer.

7. Severability Clause

Should any or more provisions of this agreement be or become fully or partly ineffective or unenforceable, it shall not affect the validity of the remaining terms of this Agreement.

8. CE Conformity

The product meets the requirements of the European Directive for Medical Products 93/42/EEC. Based on the classification criteria for medical products according to Appendix IX of the guideline, the product was classified as Class I. The declaration of conformity was therefore created by Ottobock with sole responsibility according to appendix VII of the guidelines.

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